

Product datasheet for **SC402168**

MIR6732 Human MicroRNA Expression Plasmid (MI0022577)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVIR)
Precursor:	AGGCCUAGGGGUGGCAGGCUGGCCAUCAGUGGGCUAACCCUGUCCUCUCCUCCAG
Mature Sequence:	UAGGGGUGGCAGGCUGGCC
Full insert:	>MI0022577 GTTGCAAGTGGCTGGCCCTCCCCACCGTGGGAGCAGTGCTGCCAGCTTCTGGGCGCTTCCTCACTT TCAGGAGATGGAGCCTCAGGGAAGTCCAGCGGGGCTGGCCATCCTGGAGCTGTGGGAAGTGGCACTG GGAATGGAGGCCTTGGCCTTGGGCCCCACCTGCAGAGGGCCAGGCCTGAGCCCTGGGGAAGGGGAAGG GGCTGGTGAATCCAAAGGGGAAGGGAAGTGTCTACTACCAACCAAGCCAGGGCAGGGGTGGGGTGTGA CCACCTCCCTCCTCTGCTCAGGCCTAGGGGTGGCAGGCTGGCCATCAGTGTTGGCTAACCTGTCTCT TCCCTCCCAGCCATGGGAACAAGGAGGTCTTCTCTGCCGGGGCATCTGCTGGCAGTGAAGTGGTTTCT GGAGCGGGGCCACACAGACATCACAGTGTGTGCCATCCTGGAGGAAGGAGCAGCCTCGGCCGACGTG CCCATCACAGGTGAGTGGTGCCTCTGGAGGTGGGATGGTCTTACTGCCCCAGCAGCCCTGGTTCTCCCA AGAGAGACCCTCTGGGAGTGACCTACAACCTTAGAACTCAAGCAGGGACCAGCATTTCTTTCTCAGTT TTTTCTGTCCTTAGCAACCTGCTCTGAATG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0022577



[View online »](#)