

Product datasheet for **SC402188**

MIR6752 Human MicroRNA Expression Plasmid (MI0022597)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	AUGGAGGGGGUGUGGAGCCAGGGGGCCAGGUCUACAGCUUCUCCCCGCUCCUCCCCCAUACUCCAG
Mature Sequence:	GGGGGUGUGGAGCCAGGGGGC
Full insert:	>MI0022597 GCCAGCCAGCCAGATGTGGTCAAGTCTGCTGCTGACGGACGCAGCTGTGGTGTCCCCGAGCCCCGCT GTGATATGCCCCATGCCCTGCAGGTGGAGAGCCCTGGCACGTACCAGCAGGACCCATGGGCCATGACAGA CGAAGAGAAGGCAAAGGCAGTGCCACTTATCCACCAGGAGGGCAACCGGTTGTACCGCGAGGGGCATGTG AAGGAGGTGCTGCCAAGTACTACGATGCCATTGCCTGCCTCAAGAACCTGCAGATGAAGGTACTGCCTG GAGGCTGAGGGGGAGGATGGATGGAGGGGGGTGTGGAGCCAGGGGGCCAGGTCTACAGCTTCTCCCCGC TCCCTGCCCCCACTCCAGGAACAGCCTGGGTCCCCTGAATGGATCCAGCTGGACCAGCAGATCACGC CGCTGCTGCTCAACTACTGCCAGTGCAAGCTGGTGGTCGAGGAGTACTACGAGGTGCTGGACCACTGCTC TTCCATCTCAACAAGTACGACGGTGAGCACCGGGCCCTGGGCTGCCGGGGGCTGCGAGTGGTCAGAGAG TGGCCTTTCTCTGTCACTGCTGGGGTCAAGACCTAGCCTTTCAACCCCCATTCTGAGCTCCCACGGG GGCCTGACTAAATGCCTCTACTCGGCAGGGCTGTGGGCCCC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0022597



[View online »](#)