

## Product datasheet for **SC402263**

### MIR6822 Human MicroRNA Expression Plasmid (MI0022667)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | MicroRNA Expression Plasmids   |
| Species:                  | Human  |
| Vector:                   | pCMV-MIR (PCMVIR)  |
| Precursor:                | UGGCCAGGGAACCCAGUUUGGGCUUCCGCUCUGCAGAGGCUCUAAACUGGCUUCCUGCAG   |
| Mature Sequence:          | CAGGGAACCCAGUUUGGGCUU  |
| Full insert:              | >MI0022667 GTGCACAGGGTCTGTATGTCCAGAGGAGGGAACACTAAAGGACAGAAAGAGATGGGAAAGCTTCATG<br>GGTTCCAAGCTGGAGGGACTTTTGCAGGTCAGAGCAGGCCAGGAATGGGAGTGGCGGAACAGGAGGGCTA<br>CAGGGAAGGGGTTACGCCAGTGTCTGGGGTAGCTGAGAGGAGGGCACACAGTTGGCTTCTTGC AAAGGT<br>TGGGGTACCCCTGGGGTCCCTTGGCAGTTTGGTTTATGGCAGCTCACTTCTGTCCCTTGCTGCTTCTGG<br>TTCTGGTCCCAGTCCCACCTGGCCAGGGAACCCAGTTGGGGCTCCGCTCTGCAGAGGCTCTAACTGGC<br>TTTCCCTGCAGTCTGTACAAGGCCTATGAGTACATGGGCTTCATCATGGAGAAGGAGCAGTCTACAA<br>GGATGCAGTACCAACTACAACTGGCCTGGAAGTACAGTCATCACGCCAACCTGCCATTGGTAAGGCA<br>ACCAGCCAGGGTGCACGTGAATGGACGTGGGAGGGAGGTGGGCAGGAGGGCCCCACCATGACCCAGAA<br>CTCAAGCCTGTACCCTGGCTGTTATGAGAACCAGAGGGGTTCTCAGGTCACCAGGGTGCCACCATCTTT<br>GCTCTCCATGTCCCGGTAGGCTTCAAACCT |
| OTI Disclaimer:           | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.  |
| E. coli Selection:        | Kanamycin (25 ug/mL)   |
| Mammalian Cell Selection: | Neomycin   |
| Components:               | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube  |
| Also Contains:            | Both mature and minor (star) miRNA can be over expressed from the construct.   |
| RefSeq:                   | <a href="#">MI0022667</a>  |



[View online »](#)