

Product datasheet for **SC402364**

MIR718 Human MicroRNA Expression Plasmid (MI0012489)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Human
Vector:	pCMV-MIR (PCMV MIR)
Precursor:	GGCCGCGGCGCAAGAUGGCGGGCCGGGCACCGCCCUUCCGCCCGCCGGGCGUCGACAGGC
Mature Sequence:	CUUCCGCCCGCCGGGCGUCG
Full insert:	>MI0012489 CTGCCGTGTGCAGGGGGCTAGACGTGGGGCGGGCAGGCAGGGCTCCTGGTAACAGCCCTGCAGGCCGAG TGGAGAGCAGGGTTCGGCAGGGCCGCCAGGAGCTTTCGGAAGGCCCGCCCGGCCCTTTCGAGCA GCCCGGGCCTCCGCCCTGCCCTGTGCCCAACGCCGGGAGCCCGCTTCGTCTCCAGAGCCCGCCCG GGCAGCCCGGGAGGCCGATCGCCGCTCGCGGAACCCGCGGGACCCGGGCCCTCCCGCGCGGGGCGC CCCCGTGTGACCCAGCGCGCGGCCGCGCGCAAGATGGCGGGCCCGGGCACC GCCCTTCCGCC CGCCGGGCGTGCACGAGGCCGGCTCGAAGGGGAAGTGAGTCAGTGTCCGCGGACCCGGCCGGCCAGGC CCGCGCCCGCCGCGGCCCTGAGAGGCCCGGCAGTCCCGGCCGGCGGCAGCCATGGCCGGGGGC CGGGCCCGGGGAGCCCGCAGCCCCGCGCCAGCACTTCTTGTACGAGGTGCCGCCCTGGGTCATGTG CCGCTTCTACAAAGTGATGGACGCCCTGGAGCCCGCCGACTGGTGCCAGTTCGGTGGGTGGCGGGGCT GCCGGGGGCGGGAGGCCGCGGGCTCCTGGCGCCGACGC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0012489
Synonyms:	hsa-mir-718
Locus ID:	100313781



[View online »](#)