

Product datasheet for **SC402368**

MIR762 Human MicroRNA Expression Plasmid (MI0003892)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Human
Vector: pCMV-MIR (PCMVMIR)
Precursor: GGCCCCGGCUCGGGUCUCGGCCCCGUACAGUCCGGCCGCAUGCUGGCGGGGCGGGGCCGGGGCCGAGCCCCGCGCGGGGCC
Mature Sequence: GGGGCUGGGGCCGGGGCCGAGC
Full insert: >MI0003892 GGTGTCAAAGCTTTAGGGGGCGGGGCACAAGAGGGGGCTCCCGCTCTGGGAGGACACGGCTGGGGATCAC
 GCTTTGTAGACCCCAGGAGTAGCCAGGCGAACGGGGACAGAGCCCTGTGGATCCAGGGCAAAGCTAGTG
 GGTGGAGATACACCGAACACCCGGGGCCAGAGCTGGCGAGGGCAGCGTAGACGCCTTTGGTGTGGTGG
 TCCCGCTGTGTCCCGTCCCTGGCCCCGAGCAGCGCTACCATCTCCGGACCTTCTCGATGGTCGCCATC
 ACCTTCTTGATGTATCCTTGGCCCCGGCTCCGGGTCTCGGCCCGTACAGTCCGGCCGGCCATGCTGGCGG
 GGCTGGGGCCGGGGCCGAGCCCGCGGGGGCCGCCTCCCGTCCGGCGGGCTCAGGCTCCGCGCCAGGCC
 CGGGCGCCGCGCTCTCGGCCTGCCGTCCCCCTGGAGTTGGCCGCGCCCTCTCTCGTCCCCAACGCCTCC
 GCAAGTCCCCGCCCGTTTCTGCCGCCGCGCCCCGCCCGGGCTGTCCCCGAGCGGGCCGGTGGCGCG
 GCGGCCAGAGCCAGGCCTGGAAGCCCGCGCGCCCCGCCCGACCCGGGCCCGCCCCCTTGAGCGCC
 TGGACGCTCAGCGCCCCCTGCCCGCCGCGAGAGCGGGAGTTGGGGACGGACG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0003892](#)

Synonyms: hsa-mir-762

Locus ID: 100313837



[View online »](#)