

Product datasheet for **SC402443**

MIR1247 Mouse MicroRNA Expression Plasmid (MI0004120)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CGGCCCGGGGCGCACCCGUCCGUUCGUCGCCCGGACGUUGCUCUCUGCCCCGGGAACGUCGAGACUGGAGCGCCCGAACUG
Mature Sequence:	ACCCGUCCGUUCGUCGCCGGA
Full insert:	>MI0004120 CCGCCGACGCCGCGGTTCCCGGTCAGCTGCCCTCGGCTCCAGCCGAGGTGCCGAGGGCAGGGCCC GCGGGTGGGCTGGCGCCAGCTCGGGAGACCCGAGAGTCGCGGCAAGCGACGCCGGGCCTCGCCCTCGG CCGGCTCCGGGACGCCCTGGGCGCGGGGCGCTGCGGGCGCGGGGCGCGGGGACGCGGGCGCGGGGCGGGC CCGGCCCCGGGCGGCTCTCCACTAGGGGTAGCTGTTGCCTGAACGCGGGAGCGCTCCCCCTCGCCGCT TGCCTCGCCCGGCGCAGCCCCGGCCGCGGGGCGCACCCGTCCCGTTTCGTCGCCGACGTTGCTCTCTGCC CCGGGAACGTCGAGACTGGAGCGCCGAAGTGAAGTGCCTACTTTCAACAGGGAACGTGTGCTGAGGCTTTTAGGG GTGCGGGGCAACTCCCCAACAGCGTGGGGCTCTGGTTCGGCTCTGCCACAGCGGCGCAGACCTCCTAATCT GTCTCTGCGCGTGGCGCACTGCTGGAGACTTGGGTTTCAGGGGCGCGGGGCTGCTTGACCTGGGGTCTC GAGACCCTAAGCAAGTTCCACTTTAGATACACCTCACAAATGACCCCTCAGGT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0004120
Synonyms:	mmu-mir-1247
Locus ID:	100526515



[View online »](#)