

Product datasheet for **SC402453**

MIR133c Mouse MicroRNA Expression Plasmid (MI0021863)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVIR)
Precursor: AGGCUCAUGAAGACACCAAAGGAGACAUCCCCAGCAGUCGGCUGGCUUUUCAGUUUGGUCCCCUUAAGGAGUCAGGG
Mature Sequence: UUUGGUCCCCUUAAGGAGUCAG
Full insert: >MI0021863 GCGAGAAGCCAGCAGCTTGTTCCAGTCTGACTCAGTTAAGCTCTTACAATATTTGACTGTTGGAGGAAA
ATTAATAATGGGATGTTAGGCAGTAGATTAGCTGAGACTAAATATCAACACAGCTAGGTTTATTTCTGGA
TGACTGAGATCATCTTCCGAGAACAGCAGCAGCCCTCATTATACCACAGCTACCAACTCACACAGCAG
CAGGACCACTCCCAGCAAGCAGAGCATGCATGCTGCACACAAGAGGCACAGTAGCTACAGTTTGGTGTAC
AGGCCACAGGTGATGACCTAAGGCTCATGAAGACACCAAAGGAGACATCCCCAGCAGTCGGCTGGCTGGT
TTCAGTTTGGTCCCCTCAAGGAGTCAGGGTGGAGGTAGAGGGGTAGAGCTCTACTCTGCCAGTAGACCA
AGCATATCAGTGAGGCCGACCACAATTTCTGTCCTTTCTGTATTTACTAGGCACCTACTTGCCTTAGCT
CATCACAGCTTGATGCCTGTGCTACCCTGGTGGTAAAGAAGGTCCTAGTGTGTCCTCCCAGGTCAAATG
TCCTGGGTGGATGGTGTGGCTCCTACCTCTGAGGGGCTTGAAGCCCTGATGCTTACAAACCAACAGAGT
TTCTTTGGCTTTGGATGAGCCAAGCAGACTCCAAGTCTTAAGATCAAAG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)
Mammalian Cell Selection: Neomycin
Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq: [MI0021863](#)



[View online »](#)