

Product datasheet for **SC402459**

MIR1906-2 Mouse MicroRNA Expression Plasmid (MI0015984)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GGACAUUAGGAGCAACCUCCUAGGGUUGUUGAGAAUUAUAUGAACUGCAGCAGCCUGAGGCAGGGCUGGGCAGAGACC
Mature Sequence:	UGCAGCAGCCUGAGGCAGGGCU
Full insert:	>MI0015984 TTCCTGTACATCTGTGCTCGTGTTCATCTGCTAGTGAAGTGGAGTGCTGCCCTCCCGAGGAGGGTCGTC CCTTGTGACTGATCATGCTGTCCTAACAAATGCTCTGAGCAAAAGGGTCCCTTTGGGAACCTCTCAGGAGG GGGACCCGGGTGAGGGGCGACCAGCATCTTGCTGGCAACTCCGTGGGTGGGGTGGGGTGGGGTGGTTCCT TCTGGAATGAGCACGTGGCTGACCCCCCAAGGCATGTCCCCTCCCCTCCTCCACCCACCTTCTCGGAGA TGTCCTCTTTGGGGTAGTGGGGACATTAGGAGCAACCTCCTAGGGTTGTTGTGAGAATTAATGAAGTGC AGCAGCCTGAGGCAGGGCTGGGCAGAGACCTCAAAAAAAAAAAAAAAAAAGAAAAATAATCAAAT ATATATTCTATTTGATTAAAGGCCACTCCATGATGCCAAGAGCATTTGACTAAGTAGATTATAGTGGTT AGAGGAAAATTCTTTTACTATTAATTTAATAAATGAGTACAGTAATAAAAAGAACCTTAATAATGTATTT CTATACCATGGATCAGTGTAATGTTCAAATTCATCAGAGAACTTCTTAAATAGATTGTTATTAACA CGGAGACTTACAAATGGGTAACATAGATCAGACAGTGGGAGATTTTGGAG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0015984](#)

Synonyms: Mir1906 Mirn1906 mmu-mir-1906 mmu-mir-1906-1

Locus ID: 100316809



[View online »](#)