

Product datasheet for **SC402472**

MIR28b Mouse MicroRNA Expression Plasmid (MI0016979)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVMIR)
Precursor: GACUAUAGAAAAGGAGCUCACAAUCUAUUUAGUGGCAGGCGCACACUGCCAGAAUGUGUGAGGCAUCUUAUAGUGUCU
Mature Sequence: AGGAGCUCACAAUCUAUUUAG
Full insert: >MI0016979 GGAGAGAAAATGAAAATGAAAACAATAGCATTTTTAAAGGTAGACTGTCTAATCTTTGCTGCAGTGCAAT

GGTTTTTAAAACAGAGAAATGAACCACATCGACCCCACTATGCCCATGCTGAGTATTATACAAAATGTT
 GGCAATACTGAGATCAAGAATCACTCATTTGGAAAGAATAGCAAATGCCACATACAGCACTCCTTGGCA
 CCAGCCACTGTGCTTGACTGTGCTGTACACCCCTGCTCTTTAAGTTAAGATCAGCCAAACAAGATCTTG
 CTGGCTCTTCTGGTCACATGGACTATAGAAAAGGAGCTCAAACTCTATTTAGTGGCAGGCGCACACTGCC
 AGAATGTGTGAGGCATCTTAATAGTGTGCTGATGACACTGTTTTGAAAGGAAGGTTTGAACACTTGGA
 GGTGGGGGTAAAAACAACCTAGGGAGTCTAGAAAAGCATTTTACACAGAATATTCCAAGCAAAGTAGGA
 AAGGACAACAAGAGTTCATTGGAGGCAGGATGTTGTTTAATTGCACATATTGTATAGAAAATAGTGCAGG
 TATGTACAGATGTATAGAGATGTATAGAGATGTATAGATGTATAGAGAGGGATGAATGTTGATGCTTACT
 TCACTGTCTCTTTGTTTAGTCCAGGACCACAGCTCATATTCAGATTGGCT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0016979](#)

Synonyms: mmu-mir-28b

Locus ID: 100628574


[View online »](#)