

Product datasheet for SC402482

MIR3062 Mouse MicroRNA Expression Plasmid (MI0014024)

Product data:

Product	MicroRNA Expression Plasmids
Type:	
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GCCAUGGUGAGGAAGCCAGGGAGAAUGUAGUGUUACCGUGAGCCUGUCAGCAAAUACCACUCACGGGGACACACUUUCUCCUUAUGGGUCUCCCCACUGCAGGC
Mature Sequence:	GGAGAAUGUAGUGUUACCGUGA
Full insert:	>MI0014024 CTGACCCTATCTCCTGAGTCTTGCTGAATCCTCCATTTCTATCTTCATGTTCTGAACTTTGGCCGTCTA CGACCTCTCATAGTCCCTGGTGCATCTTCCATGAGCTCCTCCTCCAGAACCCATCTCACCCGAGTCATCC GGACTCACCCCTTTTCCAGCATTGCCTGTGCTGTTTGCGGAAGAGCCCGGTCTGGTGTGGAGGTACAG GAGGCAGATGTGGCCGAGTACGGCAGAGGTACGAGAGTGTGGCCTACGGTGTCTGGAAGTGGGCCACA CAGGCGAGGCTGGACCCAGGCCATGGTGAGGAAGCCAGGGAGAATGTAGTGTACCGTGAGCCTGTCAG CAAATACCACTCACGGGGACACACTTTCTCCTTCATGGGTCTCCCCACTGCAGGCCCGAATATCAGTGAA TAAGGCTGTGGTAGTAGAGGAGCCTGTTGGGGAGTTAAGAGCCCTCTGGGAGGAAACAAGTTTCCAGCTG GACCTACTGCAGGCAGAGCCACGCTGTGTGATAGAGGAGAAACAGGGCCTGAAGGAGCGGACAGGGCCCCA GCTACTACCTGCCGCCACCTTTCTGTGGCCTCTGTGCCCTGTAAACCTGGTAAGGGAGCGTGACAGAG GGCTGGGGTACTTGGATCTGCCCTTGGGAAGGCATCTGGGGTTTTAGAGTGGATAGGTAGCATTCTGCTC AGAGA
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase .
Disclaimer:	The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0014024
Synonyms:	mmu-mir-3062


[View online »](#)

Locus ID: 100526545