

Product datasheet for **SC402487**

MIR3067 Mouse MicroRNA Expression Plasmid (MI0014029)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVMIR)
Precursor: GCCCCAGCUUCCCCAGUUCUCAGGCCCGCUGUGUGAAGUUAGGUUACCAAGCGGCUGCCUGGGAGAGGGGAGGUGUGC
Mature Sequence: AGUUCUCAGGCCCGCUGUGUGU

Full insert: >MI0014029 CCCCTTGCTTCTGTGGCAGGCACTTCACTTACAGAAACAGCCCATCTTGCCAGCCCACTGTTCTGTCAT
 GCCAGGATCTTCCTTACAAGACACCAGACTTTTCAGAGTGGAGAGTAAAGCGACAAACCTATATTTTCGTC
 TTCTCCTTAGTCCATGCTACACTGCATTTTTCTAGATTTGGGGGCTTCATTGTCTTTATGATTTAACTCT
 CCACAGGTGGTCATAGACACCGAAGGATTTCCAGAAATCCATCCACAGGGCTCAAGCCCCCGGTGCGCG
 TCTCTTTCAGCCCACTCCTGCCCGAGCTTCCCAGTTCTCAGGCCCGCTGTGGTGTGAAGTTAGGTTCA
 CCAAGCGGCTGCCCTGGGAGAGGGGAGGTGTGCCTCCTCTTCCCCCTGCATCTGCAGAAGTCTGATTTA
 TGGCTGGAGGGGTTGGCCTAGAGCGTCTTAATTTAAAGCCCCTGTCAGAGGTCGTTACACAAGCGTGGTG
 CCAGTAAAAGCTGCTTCCCTCTCCCTTTGTTTCTTTGTGAGCTCTTTCATGGCATGAAATATGAGGGAGG
 AATGAGGCTACCAAAGTGGATCTGAGAGATCTACCCCAAGCCCCCGGTGCAATCCCCTTGAAGCCGTGTT
 GGTGTTTGCCTGCCAAGATCTAAGGAAGGTGGAGGGGTTAACTTTAGGGGCCAT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0014029](#)

Synonyms: mmu-mir-3067

Locus ID: 100526532


[View online »](#)