

Product datasheet for **SC402491**

MIR3070b Mouse MicroRNA Expression Plasmid (MI0014033)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GUGCUGAGUGAGGCGUGAGCCCCUGACCUUGAACCUUGGAUCCUCCCCAUGAUCUUCUGGUGCUAUGGUCAGGGGUAGAUUCCUUGUCAU
Mature Sequence:	AGCCCCUGACCUUGAACCUUGGA
Full insert:	>MI0014033 TTAAGAAAGGGCCCTTGGGCGTGAGGGCCCTGCTCATTGCACCTGCCCACTATCTTTGAGGTCACAGTGG GAGGGGGCCCTTGTTTTTATACTCCTAGCTTGTGTGCAGATTTCTGGGTAGGAAGTATAGAGTTTGTGA GTAAGGGACTATTTGCTCAATTTCCCTCTACCCCTAGGAGCTCCTTCCCATCCCAGAGTATCCCTCAAG ATAGCTTTTTTTAAGGCAAGGATCTCAGGGGCCAACTGTGGGAATGGTTCCAGTGGGGGGCCACTGGAGC CTCCTTAGCAACAGGAGGAAGTGCTGAGTGAGGCTGAGCCCTGACCTTGAACCTGGGATCCTCCCCATG ATCTTCTGGTGCTATGGTCAGGGGTAGATTCTTGTGATGAGGCCCTTAGGACCCATAAGGGCACAGTG ACCCAGTGCAAGGCTGGAGGGGGCGTGGAGGGATGCAGTCTTGGGTCTAAGAACTGTTTTCTGGGATC CAAGCTCAAACATCCTGCCTGTTGACCCTCCTCAGCTGTGGGCATTGAGCACAGTCTGGCCATTCTCAT TTCTCATGTGACTAGAGGGCTGGGACCAGATGAGGGCTGGGAATTCTGGGTACCCTCCTTCTGCCCCATG AAACCCTTGGTCCATTTCCCAGAGGAGCCCAAAGTTTGAGAGTGAATTGACTTCCCTC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0014033
Synonyms:	mmu-mir-3070b
Locus ID:	100526519


[View online »](#)