

Product datasheet for **SC402520**

MIR3097 Mouse MicroRNA Expression Plasmid (MI0014090)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMV MIR)
Precursor: ACACAGGUGGGAAGUGUGUCCAGUGGUGUGUCUACAGAUCUGCUCAGACCUUUCUACCUGUCAG
Mature Sequence: CACAGGUGGGAAGUGUGUCCA

Full insert: >MI0014090 TTCCCCAAAGTAAACAGCTGTCTCACAGAGTGGCCTCATGCCTGTGTCTTGTGCTCTGAGCATTGTGCACTTTTGCCTCCCCTGTCTGTCACACTTTGCTGGCAGATCCTCCTCTGAGGCCCTTACAAAGCAGCCATTAGAAATTTGCAGCCATTTAGAAATTTGTAGAGGAGCCGATTCAGGAGTTTCTCCTGGCTCTTTCTTTTATATTTACTTTTCATAAGCCAGACTATAGACAGGGCCTGTGACTCAGCTCTCCAGCCCCGAAGCAGAGTCA GCCTTAGGAATGTTTGAGTGACACAGGTGGGAAGTGTGTGTCCAGTGGGTGTGTCTACAGATCTGCTCAG ACCTTTCTACCTGTCTCAGCTCTCTTACAGGGCCCGGAGATTTTCAGTAAATTCCTGTGCAGCAAAGCCAC CACTCCTGTCAACATTGACAGCCAGGCCAGCTGGCAGATGACATTCTCAATGCACCCCACCTGACATG TTCAAGGAACAACAGCTCCAGGTGACTTTGCAGACAGTCTTCCATGGGGCTTCCAGCCAAGGCTTGGGGG TGCAGAGTCACTCAGAGTAGCACGTGCCCTTCTGGGAAGGGGGCTGCTGGTCATGTGGGCATGGTGACAA GCCCATCTGACCACTGGGCTCCAGTCTCAAAGGATT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)
Mammalian Cell Selection: Neomycin
Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq: [MI0014090](#)
Synonyms: mmu-mir-3097
Locus ID: 100526557



[View online »](#)