

Product datasheet for SC402526

MIR3102 Mouse MicroRNA Expression Plasmid (MI0014099)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GUGAGUGGCCAGGGUGGGGUGGGUGGUCAGGCAGGAGAGCCAUUGAUCUAUUGUGGGUCUACUCCUGCCCCAGCCAGAGCACCCCAUUGGCUACCCACAG
Mature Sequence:	GUGAGUGGCCAGGGUGGGGUCUG
Full insert:	>MI0014099 CGCCTGCTGCCCTCTGGTTCTGCTCCCCAAGCACCACACAATGGGAGCCCTGCATGAATCCCGATTCCA CATAAGTGGGGGAGAGCCCTCACCCCTCCCTTGATTGCTGTTCTGGGGAAGTGGGTAATGTAGGTC GGCCTCTAGACTGTGCACACTGAGCCTGAATGCCACTTTGGACCCTACAGAAGGCAATTGGAGGCAAGAA GGACCGATCCCTCTTCTGTTACAGACCTCATTGTCTGTACCACCCTGAAACGAAAGTCGGGCTCGCTG AGGCGCAGCTCCATGAGCCTGTGAGTGCCAGGGTGGGGCTGGGTGGTGCAGGCAGGAGAGCCATTGATC TATTGTGGGCTCTACTCCCTGCCCCAGCCAGAGCACCCATTGGCTACCCACAGGTACACAGCAGCCAGC GTCATTGATACAGCCAGCAAATACAAGATGCTGTGGAAGCTGCCACTGGAAGACACAGACATCATCAAGG GTGAGGCTGGCTGAGTCTGTCAGGGTTCTGCTGGTGGCTTTCCCTGAGTCTGCCTGCCTGTCTGTTGG TTTGTCTGCTCTTGTGACTCATCTTTCATCCCATATGAGAGATTGGTTTTGAAGAACCAGGGTTAGTTTCC ATGCAAATCATCCTCTTAGTGCCCTAGCAAGTTCTCCAAGAGTTGGCCGTTCTGTGCCCTGTGACTCTAA ACCA
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase .
Disclaimer:	The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0014099
Synonyms:	mmu-mir-3102



[View online »](#)

Locus ID: 100526508