

Product datasheet for **SC402548**

MIR3470a Mouse MicroRNA Expression Plasmid (MI0014696)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMV MIR)
Precursor:	CCUGGACUCACUUUGUAGACCAGGCUGGCCUUAAAUUCACAGAGAUCAGCCUGCCUCUGCCUCCUGAGUGCUGGGAU
Mature Sequence:	UCACUUUGUAGACCAGGCUGG
Full insert:	>MI0014696 AGATTCTAATTAATCAGATGACAGGCTCACACCATCACACCTTCTTCTGGTTGCTGGTTTGTATCCA AGGCTGGAGTGT AAGGGGAGGTGGTACCAGCTGGACAAGAAAGT AAGGAACACTGTTGATGTTGATGGCT CAGTACCAGCTAACAAGACTCCGTTGCCTCAGCTAAAATGATTTGTCCTGAAGCTGTTGGGGCATGGGC TGGGGAGAGGTGGCAGTTTAGTTTGAGGGCATTATTTGTTTATGGTTTTATTGAGACAGGGTTTCTTTG TGTAATAGAGCCCTGGCTGTCCTGGACTCACTTTGTAGACCAGGCTGGCCTTAAATTCACAGAGATCAGC CTGCCTCTGCCTCCTGAGTGTGGGATTAAGGCATGGTAACCACACCAGTAGCTCAGGGATTATTCTT CCATTAGTAACCTTAAGAAGTGTTCATCCAGCATTCTTGGTTTGTGGTGGTTGGAAGGTCAAGCTGCG TGCTTCAGGGGCCCTGCTAGAAGTCTCCTGGTCTCTCAAGAGTGCCTCCTTTCACTGTGGGGTATTTTA GCCTTCGGCTATGTCCAGCTTGTCTTACCCTGTTCTGGGGACTAAACCCTGTCCATTTAGCTAAATC TCCAGCTCCATGGACAACCTTTTACGTTATGAAATGGCACTGCTGT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0014696
Synonyms:	mmu-mir-3470a
Locus ID:	100499518



[View online »](#)