

Product datasheet for SC402556

MIR3473d Mouse MicroRNA Expression Plasmid (MI0018033)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVMIR)
Precursor: UGUGCUUAACCACUGAGCCACUUUCCAGCCCUUUGAAAACUUUUUGAGGGCUGGAGAGGUGGCUCAGUGGCUAAGAGCAC
Mature Sequence: CCACUGAGCCACUUUCCAGCCCUU

Full insert: >MI0018033 TTTCTGACTATAGGTAAAAGTGCTTGAGAATCTCTGTGTGACTCTGCTACCCAATCCAAGGTCTCCTTGTG
 TCAAGGAAATTTAAAATGTCATTTTAAGCTTTATTCACCTTTATTTTAAAAAAAATTTAAATTTCTGT
 GAAAGCGTTTTGCCTGCATAGGTGTGTGTGTCTCTTGTCTTACCTATGTTCAAGGATGTCAGAAGAAAG
 TTTCAGATCTCCTGGAAGTGGAGTTAGACAGTTGTGAGCCACCATGTCAGTGTGAGAATCAAATTCAGA
 ACCTCGGCAAGAGAAGCAAGTGTGCTTAACCACTGAGCCACTTTCCAGCCCTTTGAAAACCTATTTGAGG
 GCTGGAGAGGTGGCTCAGTGGCTAAGAGCACTGACTGTTCTTCCAGAGGCCCTAAGTTCAATTCACGCA
 ACCACATGGTGGTTCACAAGCATCTGTAGTGGGCTCAGATGCTGTCTGCTGGTGTGCATGAAGACAGAGC
 ACTCATATACATAATGTAAATAAATAAGTTTTTTTTTAGAAAAAAAATAGAAGCCAGGCGTATGCCTTT
 GATCCCAGCACTCTGAGTTCGAGGCCAGCCTGGTCTACAGAGTGAGTTCCAGGACAGCCAGGGCCACACA
 GAGAAACCCTGTCTTGAAAGTCAAAGAAAAAAGAAAAAAGAAAAAACA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0018033](#)

Synonyms: mmu-mir-3473d

Locus ID: 100628592


[View online »](#)