

Product datasheet for SC402580

MIR3968 Mouse MicroRNA Expression Plasmid (MI0016977)

Product data:

Product	MicroRNA Expression Plasmids
Type:	
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	ACGCCGAGCGUGUGGUGGUAGGAUCCGUGUAGCACCGUCUGUGUGUUGUGGCCACACAGAGAACUUGAGCGAAUCCACUCCAGACACCAUGAGUAGG
Mature Sequence:	CGAAUCCACUCCAGACACCA
Full insert:	>MI0016977 TCTGCCCTGGGAACAATGACCCTGGCACCTTGTGCCTCTATCCCCAACCCAGATTACTGCAAACCTGAA CTTTGACCCCTGCACAGCAGTACATAGAGTCTAACTTTATTGGCAGCAAGCCTGGGCCGGCGGGCGCTAC CTATGCCTCAGGATACTGCTTGAGTGTGAGCGAGGCTCCACTGCTGACCATGACCGCGGGGTACAGGGGT TCGAGGAAGGAGGTGAGGAAGCGGTAGAGGAGGTTACGTCGCCGCCACGCTGTAGAAGGACAGGCAGC CAGTCTGAAGTCCAGCGCCACGCCGAGCGTGTGGTGGTAGGATCCGTGTAGCACCGTCTGTGTGTTGTT GTGCCACACAGAGAACTTGAGCGAATCCCACTCCAGACACCATGAGTAGGGGTTGCGGCCAGCAGGTAG ACGATGTTGCGACGGGGGCGACCGCTCAGCGGAGGGTGCAGCGGTAGGTGACTCCCAGGCACATCCACG AGTTGGACACCTCAGCTCCAGTAGTGACAGCCTTCGGACAGGCCTGGGCAGCACAGAACCTGGGGCCA CTGCTTGAAGCCCGTAACGGGCTGTTGGTGGTCGAGGGCTCCAGAAGGATGCCCAAGTTCAACACCGAG TGTGTTTCCTGAACAAGAACAGCTCTTCGCTGGCCGTGCAGGGGTCAAAGTTCAGGTTGCAGTAATCTG
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in
Disclaimer:	miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0016977
Synonyms:	mmu-mir-3968 mu-mir-3968
Locus ID:	100628623


[View online »](#)