

Product datasheet for **SC402595**

MIR466o Mouse MicroRNA Expression Plasmid (MI0014052)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVMIR)
Precursor: UGUGUGUAUGUGUUGAUGUGUGUGUACAUGUACAUGUGUGAAUAUGAUUACAUAUGACACACUAAGACACAUUGAGCACA
Mature Sequence: UGAUGUGUGUGUACAUGUACA

Full insert: >MI0014052 CCCAGTGTATGTGCCTGTGTGTATATGTATGTGTGCATGTGCATGTGTGTATATGAGTATACATATACAT
 ACGCACATACCCATACAAGCATGCACACAAACCCAAAACCTGCCACTCATCGATCTATTTGAATGATATGT
 GTCTTAGAGTCTTAGAGATTCTGTATGACTGGGCCAGGAGAGGATTCACCAAGTAAAGACCCCTTG
 CCAGAGTCAAGATAATCAGGAGTTCAAGGTCATGGTAAAGGCATGATGAGTTAGAGGTCATCCTGGCATT
 CATGTGATTGCCCTATGTATGTGTGTATGTGTTGATGTGTGTGTACATGTACATGTGTGAATATGATAT
 ACATACATGCACATAAGACACATATGAGCACATAGACACACAGGAATGGCACTCACTCATCAATCTAT
 GTGACTGAGAGCTTTGGAGATTATGTGTGACTAGGTCAGGGAGAGTTATCACCCATCATATTAATAATGC
 TTGTTAGAATCAGGAGGTCAGGAGGTCAAATTCATGCTTTATTTTCATGATGACTTAGAAATCTGCCTGGG
 ATTCATGTGATGCCCCATGTATGTGTACATGCGTGTACCTGTGTGTGGATTGTGTATATGAATATACAT
 ATCTACACACGCATAGACATGTACACAGACACACAAACATGAATGCAACTCATC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0014052](#)

Synonyms: Mir466o mmu-mir-466o

Locus ID: 100526482


[View online »](#)