

Product datasheet for **SC402602**

MIR467a-5 Mouse MicroRNA Expression Plasmid (MI0014062)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVMIR)
Precursor: UGUGCCUGUGUGCGUAAGUGCCUGCAUGUAUAUUGCGUGUAUAUUUUAUGCAUAUACAUAACACACACCUACACACACAUGCACA
Mature Sequence: UAAGUGCCUGCAUGUAUAUGCG

Full insert: >MI0014062 GCATATAATTGTGTGTGCACGTACATGTGTATATATGATCATCCCTGTAAATACACACACATACAAATGT
 ATTCACACATACAGGCATGCACACTCACACACATAAGAATGACACTCATTGATCCATTTGAATGATATAT
 GTCTTAGAATTTTGTGATTATAAAGTAGAGGGACAGTGAATATTCTCTCTCACCAGGTGTAAACAATTG
 GCAGTGTGAGGAAGATCAGGAGTTCAATTTTATGCTTAACTTCATGAGCAGTTAGAGATCATCCTGGAAT
 AATGTGATGGCCCTATGTATGTGCCTGTGTGCGTAAGTGCCTGCATGTATATGCGTGTATATTTTATGC
 ATATACATACACACCTACACACATGCACACAGACATGCGGGAATGGCACTACGGATCCATTTGAA
 TGATATGTGACCTACTGCTTTGGAGGTTATGTATGACTGGGCCAGGGTGAAGTCTCACAACACGGGAA
 ACCCTTGGTAGAGTCGGAAGATCAGGAGTGCAAGTTCATGGTAAAGGCATGATGTGTTAGAGGTCATCC
 TGGCATTGATGTGATTGCCCTATGTATGTGTGTATGTGTTGATGTGTGTACATGTACATATGTGAAT
 ATGATATACATATACATACACGCACACATAAGACACATATGAGCACACAGACA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0014062](#)

Locus ID: 100526487



[View online »](#)