

Product datasheet for **SC402626**

MIR5121 Mouse MicroRNA Expression Plasmid (MI0018030)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMV MIR)
Precursor:	GGUGAGUUGGUGAUGUGGUGACAUGUAGGACAGGGUAGGCUCUGGCAGCUUGUGAUGAGACAUCUCCACUCAU
Mature Sequence:	AGCUUGUGAUGAGACAUCUC
Full insert:	>MI0018030 CTCCCTGTTAACTCTAGGCCGGAAGGTGGTGGTTCGTACGCTGTGAAGGCATCAACATTTCTGAAACTTC TACAGAAACAAGTGTGAGTTAGGCTGGGCGAGTGGTCTTCAGGCGGGCTGTTGCCTTGCTCCTGGACAA GCAGGCTATCTGCTCTGCTCTTGGGTCCAAGGGTGCCTCACATACTTTTGTCTGACAGTAAAGTATCTG GCCTTTCTCCGGAAGCGGATGAATACCAACCCCTCCCGAGGCCCTACCACTTCCGAGCCCCAGCCGCA TTTTCTGGCGCACTGTGCGAGGTGAGTTGGTGTGTGGTGACATGTAGGACAGGGTAGGCTCTGGCAGCT TGTGATGAGACATCTCCCACTCATGTTTCGAGTTGCTCGACTATGAGATGACTCTACATGCACTACCATCT GAGGCTGTTCTAGCCCTTGCCCTGGGGTCAGACGGGGTCTCAGCACCTCTGATGAGGCCGCTGTTACC CCCACAGGCATGCTGCCCCACAAGACCAAGAGAGGCCAGGCTGCCCTGGAGCGCCTCAAGGTGTTGGATG GGATCCCTCCACCCTATGACAAGGTGAGCCTTGCAAGCCTAGCCTTACAGCCTCTGCGTCCAGCGTCTGT GATGTTCTGCTATTACCTACATTGTTTGAGCCTCATGAAAACCC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0018030
Synonyms:	mmu-mir-5121 mu-mir-5121
Locus ID:	100628629



[View online »](#)