

Product datasheet for SC402641

MIR5135 Mouse MicroRNA Expression Plasmid (MI0018047)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVIR)
Precursor:	GAGAGGACUCGGAGCCAGCCACCUAGACUUUCCAGGAGAAGCGCAGGUCUAGGUGGCAAGGGCGUCCUGGUCCUACU
Mature Sequence:	AGGUCUAGGUGGCAAGGGCGUCCU
Full insert:	>MI0018047 CCAGGAATGTTTCCAGTGGGTGGAGACTGCGCTCGCGGCTGCGATCCAGCAATTGGCGGCTCTGGGGGGC CCGGGGGAACCTTCGCGGCGGTGCGAGCGAGCTGTGGCTCCCGCTGTAGACAGCGCCGCGGTGGCTAAT GATTAACGCGCGTGGCCCGCGGTGTGGTCCCGCTGGCCGCTCCCTTGCCACCTGGAATCTGGGATC CCAAGGCGCACGTGGTGGGTTACTCTAGTGCCACCGCCTGATATCACGCCGAGAGGGCAACGAAAGCC TCCAGACCAGGTTCCGACCAGAGAGGACTCGGAGCCAGCCACCTAGACTTTCCAGGAGAAGCGCAGGT CTAGGTGGCAAGGGCGTCTGGTCTACTGGAGCAGGTACCACGACCGCTGAGTCAGAAAGGCACGTCT CGAATTCCTCGCCGAACTAGGTTTCGGTGGCTCTTTACAGTTGTAAATCTTACCTTCTGAGGAGAAAG AGATAAGAAGCGCATAACTTGGGGCGGGCGGGCGCGCGCTCAGTGGGACTGTTAATGAGAAAACTG ATCGAGGAACCTGAATGTCAGTGGACCCATTTTGAAGTGCATCTATCAACAAAAGGCATGTACCGGGAA ATTAGGGGATGGGAAAATTGGGGGATGATTAGCACAAAGCACATAGAATA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0018047
Synonyms:	mmu-mir-5135 mu-mir-5135
Locus ID:	100628604


[View online »](#)