

Product datasheet for **SC402651**

MIR5622 Mouse MicroRNA Expression Plasmid (MI0019190)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMV MIR)
Precursor:	UUCACCACACCCAGCUUAAAGAAAGUUUAAAAAAAAUUUCUUAGCUGGGUAGUGGUGGUC
Mature Sequence:	UUCACCACACCCAGCUUAAAGA
Full insert:	>MI0019190 CAGACTAGTCTTAAACTCCCATCCCTTTGACTCTACCTCCTGAGTGGTCTGCTATGCTATTACACCT GGCTCTGGTTTAGTTTTATATAAGGGTTTTGGTTTAAATGAACATAGATTTTATAAAATGAATTTTAA CAGGCATTTAGAGAAAGTTTTTTTTTTGTGTGTGTGTTTCATTTTTGAGGCATGTTTCGCTATATAGC TCTGGCTTGGTCTGGAACCTCAGTACGTAGACTATGCTGGATTCAAACTCCAAACCTGCTTTTGTCTTCTA TTGTGTGATTAAGACATGTTTCACCACACCCAGCTTAAAGAAAGTTAAAAATATTCTTAGCTGGGTA GTGGTGGTGCATGTCTTTAGTCCAGCAGTTGGAAGGCAGAGGCAGTTGGAGCTCTATAAATTTGAGGCC AACTGGGTCTACAGAGTGAGTTCCAGGACAGACCTACACAGAGAAACCCTGTCTTGAAAAACCAAACC AAAACCAAAAAACAAACCTTATTTGTGTATGTGTGCCCTGACACATACATATGAAGGTTGGAAAATTGCT TGCAGGTGTCAAATCTCTCTTTTCTTGTGGATACCAGAAATGGCACATAGGGTTGTCGGGTATGACTATC CCATGATCCCAAGAGAGGATAACGTTTATG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0019190



[View online »](#)