

Product datasheet for **SC402659**

MIR599 Mouse MicroRNA Expression Plasmid (MI0012529)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVIR)
Precursor:	UGCUGUCCACAGUGUAUUUGAUAAGAUGACAUAGGAGAGGAACUUCUUUACCUUUGUGUCAGUUUAUCAAACCAUACCUGGAUGAC
Mature Sequence:	UUGUGUCAGUUUAUCAAAC
Full insert:	>MI0012529 AATTCAAGTATGTTTATAATGTAAGTGTGCAATGTTTACCAGAATGAAAGCAAATTAAGAATTATCTCT TTCATTAATGTCTCTTCTGCCACTAAGATGAACAAGGGGTGCTTCCATCCTATAACCTGTGCATTTGTT TATTTTTCCTAGGTTTGGGAAATTTCCAGTTCTGTGGTACTATACCTCAGTTTTATCAGGTGTTTATTA AAATCACCTGAAAATACTGAGGCTATGTTTCACTGAGCATAAGCTTCATGGATTCTTGAGTGACTGCAG CTGTACTCACCACAAAGATGTGCTGTCCACAGTGTATTTGATAAGATGACATAGGAGAGGAAGTCTTTTC ACCTTTGTGTCAGTTTATCAAACCCATACCTGGATGACAAGCAATTCATGGAAGGAAAACAGTCCCAG AATCAAAGCATTATTTCTGTAGAAGAGCTCTGTGCTCTCTAAATGCCAAGGCATAAATATACGACAG TTTCACTAAATGACTAAAGCCTTAAATAATATGGATAGAAACGGTTGCACTGCTAAAGTTATAACCACAG TTATTAAGGTTTACTGGTATGCTAGGGGATTTGGGTGAAAATAGCACTAAGGATAAATATGTATTACCTG ATAATAAATATGAGGAATTTATGTTAAATATGCTTTGACTTTGCTCAAAATTATT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0012529
Synonyms:	mmu-mir-599
Locus ID:	100316736


[View online »](#)