

Product datasheet for SC402677

MIR6345 Mouse MicroRNA Expression Plasmid (MI0021873)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: GCUUUUGAGUUUUUCCUUGUUUUCAGUGUUAGUCUAGCACUGCUCAGGGUGAACUAGUACAGUGUGUGGUAGUACUAGGUUCUCCAUGGACACAGGCAUAAUAGAA

Mature Sequence: UAGCACUGCUCAGGGUGAACU

Full insert: >MI0021873 GAATCAGGGAAAGAGGAGCTTTAATCCATAACCCCTTTTAAACAGAAGTGTCTAGCGACATTCTGACT
 CAAGGACTCAGCTTTACCTAGCTATAGTCTTATCTGGACTTCACTCCTTTGTTGGTGTGGGTTTTGTT
 TTGGTTTTTTTTGGGTTTTGGTTTTTTGTTTTTGGTTGTTTTTTTTTTTTTTTTTTGAGCCAGGCCGG
 GCTGAATTTTCTGTAGCTGCAGATGATCTTGACTCTCAGGCCTCGCCTCTCCAGTCCTGACCTTATACAC
 ATGAACCACCATGTCCAGCTGCTTTGTGAGTTTTCTTGTTCAGTGTTAGTCTAGCACTGCTCAGGGTG
 AACTAGTACAGTGTGTGGTAGTACTAGGTTCTCCATGGACACAGGCATAAATGAAACTTGGTTTCTGCTC
 TTAAGCCCTAGGAGACAGTGTTCATAGTCTGCAGGTTGCCTTCACTGCCCTCCAAGCTCCAGGCAACCT
 CATGATTCCATAAGATGACATCTCCAGTGTGAGCCAGTCAGAGAACAATGGAGCCATGGTCAGAACCCAC
 TTCCTAGTGAAGATCTGGGGTCAGATACCCTCTAGCTAAGATCCCACTCCTATACAACTATGTT
 CATTCAAGACATATTTTGTCTTTCTAGGGTTCTAGATATTTCTTTAATATGCTGAGGAACATTGAAG
 GGATTGACAACTGACACAGCTGAAAAA

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect
 mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0021873](#)



[View online »](#)