

Product datasheet for **SC402681**

MIR6349 Mouse MicroRNA Expression Plasmid (MI0021877)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Mouse

Vector: pCMV-MIR (PCMVIR)

Precursor: GCUCUGAGAAUUCUUAAUCCCUCCACCCACAGAGUAAGCUAACAGCCUCCGUGUAACCUGUGGAAUGUGGGAGGGGCAUGCAGAAUGUGGGAGG

Mature Sequence: UGGAAUGUGGGAGGGGCAUGC

Full insert: >MI0021877 GGTCAGTAGGTACATGACAGCCATGTTGGCATCTCACCAAGGAAGTGTGGACTTTGGCAATAGGCTGG
AAGGAAGTGACGTTTGCAAATACAAGTGATCTTCAAACCAACTCCAGGTCTCAGGCTGGGAATGGATAG
GTCTCTTTATAGGATCCTTGAGAAAGAGGCTCAGGCCTCCTTCTGTGGCTTCCAGAGTAAAGGATCAGAA
TAGGACAGTCCAAGGAAGGGTGCTTGTTAAGCTGTAACACAGAACAGGAATATATCCTTTCTACTA
GCAAGACAGAAGCTGTGAGTGCTCTGAGAATTCTTTAATCCCTCCACCCACAGAGTAAGCTAACAGCCT
CCGCTGTAACTGTGGAATGTGGGAGGGGCATGCGGAATGTGGGAGGGGCATGTGGAATGTGGGAGGGGC
ATTCAGTCTGTGAGAAGAAACAGAAGGGAGGGAAAGTGCTGAGGGGAGGGCTGGGAAAGGCAGGTGGTG
AGAGCAGCAAAAGTGGGTGAAGTGTGAAGGGGCCAGAAGGGCAGGAGGTGGTGTAGCTCCTGTAATTAG
TGGAGAATTAGCAGAAGAATCAGGGCTCAAACACAGCCAATGGAGGAAGGATCCTTTCTGCTCAGAGAGG
CTTCTGCTCTTGCTTGGTGTGGTGAAGGCCAGGCAAGCTGTTGAGACATGCACTTCTCTAATAAG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0021877](#)



[View online »](#)