

Product datasheet for **SC402694**

MIR6362 Mouse MicroRNA Expression Plasmid (MI0021891)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: AUUCAUUUUUGCUGGUGGGCAUCUAGGCUAGUUCUGUUCUUGGGCUGUGAUGGACAGGACAACAACAAGCAGGGUGUGGAGAGUCCUUCUGUGAUGGAUAUGG

Mature Sequence: GUGGGCAUCUAGGCUAGUUCUG

Full insert: >MI0021891 ATGTCATGTGAGACCTTTTACTTTGCTTTCTCAACACTTCTTGCTCCCTCTTAAGATCTCTGATCTGAT
TTCACAGTCTATAACCAACCTCATAATCCCTCATTTTCATACACAGAACTGTAAGTTGGGATTGAGT
AGGAGAGAGAATGCATAAATTTTTTTTGTCTTTCTGTGTCTAATTCACCTCACGTAATAAAATATTTTCC
AGATTCATCTATTTTCTATAGATTTTCATAATGTCATTTTCTTTACAGCCTAGTAAAATCTATTGCTC
ATATAGGCCACATTTTCATTATTCATTCTTTTGTGCTGGTGGGCATCTAGGCTAGTTCTGTCTTTGGGCTGT
GATGGACAGGACAACAACAAGCAGGGTGTGGAGAGTCCTTCTGTGATGGATATGGAGTCCAATGGGGATA
TGTACAGGCATGCTGCAGCCAGGTGGCTCAGCCAAGCCCCACCTCTGCTCTCAGCCATGTGGTGTCCCT
CTGTGGTGTCTCAGCACTTACTAGTTCTTGGCCAGAAGAGCTCAAGCTCCAAGCCTCTCTGGCATCCCCCT
TAGTTGAAATCCTATTCCTGACTGCTGTTTCAGCTAATGCCAATGATGCCATTCCAACGCAACTTTGAT
CTCTCACTGACCTGGAAGTGGTTCTTGTCTTTGCATAATTAGTTGCAGTACCTGCCACACATACACACT
CTTCTAG

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#).

Disclaimer: The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0021891](#)



[View online »](#)