

Product datasheet for SC402703

MIR6371 Mouse MicroRNA Expression Plasmid (MI0021900)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: AGAUGGGAACUCAUUGUACAUUUUUGAAGUAGUUUGUCAUGCUGGAGGUUCAUUUAGGCAGUUCUAGUGGCUUGAAGCUCUUCUGGGAAGUUCAGUUGAUCUCUCU

Mature Sequence: UUCAAGUGGCUUGAAGCUCUUCU

Full insert: >MI0021900 AAGGAGTTCTCATGTGACTGATTTATATTACATGCAGTTTAAGTATTTCAATATGGAAATTTTGTCTTC
 GTTAATACCTACCTCTAATTCATTTTAGCTCAGAATAGGTTTAAGGAATGTACATGTTGGTTATAGGACA
 AAAAAATAATTCAACATGTGACTTTTTGTGTGTGTGAGAGTCAATAAATTTTACTGGCATTATTTAT
 AGTAATATAGGTGAGAGGTTATGGACAGGAACAGAAATGACTCAGACAGGTGCATTGCAGCCTGCTTCAG
 CTCAGGGGACAGTTCACAAAAGATGGGAACCTATTGTACATTTTGAAGTAGTTTGTATGCTGGAGGTT
 CATTTAGGCAGTTCAAGTGGCTTGAAGCTCTTCTGGGAAGTTCAGTTGATCTCTCCTATCATGTAGCTGG
 GCTGATCTGTTTTCACTCAAGTAGCTGCTCTCATCTGAACCTTTCTAAGCAACTCAGCTTATCCCAACT
 CCTAATTCTTCATGCTAGTAAGCTTGTATGGTAAATTATTTTATAATGTATTTCTTGGGCCAGTAAG
 ATGCCTTAGTGAGTAAACCACTTGCCACGTAAGCGTGATGACCTGAGCTAGGTCCCTGGACCCACAGT
 GGAAGAAGACTCACTAAGGTTGCTTTCTGATCCACACACTAGTGCTGTGCTAGGCACACAAACACACAGA
 GCATGGAAGGG

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms
 should not affect the function of the mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0021900](#)



[View online »](#)