

Product datasheet for **SC402717**

MIR6385 Mouse MicroRNA Expression Plasmid (MI0021916)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: CCAGGGCAAGGAUGACAGAGAGGAAGACUAGGGGCUUCCAGGACAAGGAGUCGCAACUCUCUGUAACCUGGGUUCUGCAGGUUCCUCAUGUAUCUGGGCCCUU

Mature Sequence: GCAAGGAUGACAGAGAGGAAG

Full insert: >MI0021916 ACCCCGCGTGTTCCTATCCGCCTTCTATTCTGTTCTACCTCCAACCCCATCAGTCCCTGCTGGATCAG
 AGCTCTGGGCTGCATCCCTTAAGAGCATTGCCAGAGACTCCTCCTGAAGCCCAAAGATAGGTTTTTTGT
 TTGTTTGTGGGTTTTGTTTTTTTTTTGTTTTTTTTTTTTGGGTTTTTTTTTTTTTTGTGTGTGTGTGC
 TATGAAATACCTGATGATTCTGGATGGCTAGCCCCAATACAGGGCCATAACTATTGGCTGAGT
 AAAGGGTCAGAGCTGAGGAACCAGGGCAAGGATGACAGAGAGGAAGACTAGGGGCTTCCAGGACAAGGAG
 TCGCAACTCTCTGTAACCTGGGTTCTGCAGGTTTCCTCATGTATCTGGGCCCTGTTCTAGGAAACGTGTG
 GTCACAGCCAAAGGTCAAGAGTAGATTACTGCCAGTCAGTGGTGGCTCACACCTTTAATCCCAGCACTTG
 GGAAGCAGAGGCAGGCGGATTTCTGAGTTCAAGGCCAGCCTGGTCTACAGAGTGAGTTCAGGACAGCCC
 CGGATACACAGAGAAACCTGTCTCTAAAAGGTCAACTATGGCCTATACACTGTTCCCAATCTCTCCACA
 GAGGTTTTGTCAGAAGTGAATCTGTGGAGTTACACATCTCCCTCTCACGTCCTGGTTCTACTCCAGGCAG
 ATATG

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#).

Disclaimer: The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0021916](#)



[View online »](#)