

Product datasheet for SC402721

MIR6389 Mouse MicroRNA Expression Plasmid (MI0021922)

Product data:

Product	MicroRNA Expression Plasmids
Type:	
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	CCUCCUAAGCAGUGCAAUGUUAACUUUGCCAAACUCCCAAGCUAGUUUUGCCUAGCUUUGGACUAGUGUGUAAUUGUGUUAACAGCUUUUGUUGGAGU
Mature Sequence:	CAGUGCAAUGUUAACUUUGC
Full insert:	>MI0021922 CTTGTGGAGAAGGATTTAAATCCAATCAGAAAGTGCTATGTTGCTGTCATGACATTTGTATCACTATTGC GCCAGTGGGGACATCTTGCTAGGCTATTCAATTATTGTAGTTCCAGTGTTATAGCTAAGGGAGATTCTTT ATGCCTTTTCTCCCTGGTGTATATAGAACTTCGAACACTATAAAAGCTAGGGATGATATTTCTGCTA AGTACCAGTAGATTTTCCCATGTCTTCCCTTTAGACCTATGAACTTATCTGTTGTTAAGTCTGGTCT CTTCCCTGGATCCACAGTTTCTCCTAAGCAGTGCAATGTTAACTTTGCCAACTCCCAAGCTAGTTT TGCCTAGCTTTGGGACTAGTGTGATTGTGTTCAACAGCTTTGTTGGAGTGTCTTTGGGTCTGTGGCC TAGCCTGGGGTCATGTTTGGCCTCATCCTTGGTTCCAGATTGACAGAAGCCACTGAGAGCTGTGGAAC TTCCCATGTGCCTGTATTATTTCTCACATTTAACCACAGGTAAGAAGGGCCACATAATATACAGGAAAAT ATGTGACATTGCAGGTATCACTACAACTTTACAGAAAATGACATGTGGCCACAATTCATGAATGACAC TTTAGGGAGTTTCTTGGGATTCTTTTGTACATAAATTTAGCTGAGTCCTATACATTAGAGAGCA GG
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in
Disclaimer:	miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0021922


[View online »](#)