

Product datasheet for **SC402735**

MIR6403 Mouse MicroRNA Expression Plasmid (MI0021939)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: AUUGUCCCAACCAGUCUGGCCUACCCACUGCCGCCUGUGCGGCCGUGCCUCUGUACCCCAUUGCAACUUGGCAAAGGCGCUGUUUCCAGUGGCAGAGCCCC

Mature GCAAAGGCGCUGUUUCCAGUG

Sequence:

Full insert: >MI0021939 CTGTCCAATGGGGAGAACAGCATGCAGCTGCCTGAGCCCTGCTATGCCCCACAAGCCAGCATCTAACTG
GCCCTGCATGAGGCTCTGACAGCAAGAGGCCAGCGACAAGTGTGAGGGCTGGAGCTGGTTCCAGGCCGGT
TATGTCCGTGCTGGGCTTGAATCGGGTGATCCATTCTTTGCACGGGGCTCCCTGCCAGGGCCATAAA
CAAACGGGAAATCTCAGGAATAATAAGAGCTTGCTGTTTGTCTGGGACAATCATCCAGGGCTGTAAAAG
ACAGCCGCTTAATCGAGCATATTGTCCAACAGTCTGGCTACCCACTGCCGCTGTGCGGCCGCTGCC
TCTGTACCCCAATTAGCAACTTGGCAAAGGCGCTGTTTCCAGTGCCAGAGCCCTGAGCCATCTGGGGG
CGGGGAGAGGACAAAGAAGCTGAGCCAGTGACTCCTAAGGAGATACCTCCAGACTGCCAACAGGCT
CTAGAGAAAATCCCCACCCCGGGACCTCTGGGGGCCAAGGCCAGGCAGGAGTGAGGACATACCTAGGA
CTAGAGTCAGTATTAGCTGGCCTGGACAGGTGGTGACAGGCAGGAGGCACAGTCTAGCTCAGCCTTGA
CTTACTGTGTGGCCTTGATAAAATGCTTAGCCCCCTCCGATCCTGTTTCTGTATGAACAAAAGATTGCCG
ATTGCCGTTTGCAAGAGAGAGAGAGGGGG

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). Th

Disclaimer: miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0021939](#)



[View online »](#)