

Product datasheet for **SC402746**

MIR6414 Mouse MicroRNA Expression Plasmid (MI0021952)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVMIR)
Precursor: UAUAAGGCUUUUCAGACCCUGCUCUUUGGGGGCUUGUAGCUCCAGGUUUUCCCCAGUCUAAGUUCUGAGCUGAGUCUUGUU
Mature Sequence: CUAAGUUCUGAGCUGAGUCUUG
Full insert: >MI0021952 TTGACTTCCTCCAATAAACCCAGCAGCAGATCTGTGCTGCCTGCGGAATGTGGCGTTAAATGCAATGGGT
 CAGCATTAAAGCTGGAACAAAATGTGATGACGAATTAGGAAATGAATAGTAAGCCATGGACGACAGATCTT
 GATTGTGATAATGGTCATTCACCTCGCTGACTGGAACAAACACAGCTGCAAGCTGTATCTTCAGATTGACA
 AAAGTGGACCCACATATCCTGGTCTTGGGCAAATTAATTCTGGCTAGAATGGTCATCTATGAGGAAATTA
 GGTGTCTGCAGTGGAATATTTATAAGGCTTTTCAGACCCTGCTCTTTGGGGGCTTGTAGCTCCAGGTTT
 TCCCCAGTCTAAGTTCTGAGCTGAGTCTTGTCTCCTGTGGGTGTATGGAAATGTCAAATGCAGGACATA
 ATGCAGAGGAAAAGGCAGCTCCCTTTTGCAGATGGAAGGCCAGAGCAAGAGCTTAAACACACTTGCA
 CAATAGATGATTAGCAGTCTGGCTAAATCAATAGGTCACGAATAGTAATTGGTGTAGAACCAAAATCTGT
 TAAGCGTTAATCTCCACCTAACCCCTCATCCCTGTATAGATAAGTTTGGTAACCTTTTCAATATCATCATC
 TCTCAAGGCCCAAAATCAAATAAACTGTCTGCACAAATGTTCTTGGATGC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0021952](#)


[View online »](#)