

Product datasheet for SC402748

MIR6416 Mouse MicroRNA Expression Plasmid (MI0021954)

Product data:

Product MicroRNA Expression Plasmids

Type:

Species: Mouse

Vector: pCMV-MIR (PCMVMIR)

Precursor: UCCUCCGUAUCAUCUGCUCUUCUGGUGGUCUGCUUACUGGGAAGACUGCUCCUCCUCUCCCCUAGAGCGUUCACCUGUGAGUCAGCCUUGGCAAAGAGCAGCAA

Mature Sequence: CUCCGUAUCAUCUGCUCUUCU

Full insert: >MI0021954 TAAGCACATCAGGGGAAAAGCACCGAGACTCCCCGGTGTGTGTTGCAAAGGTACTCGGGCAACAGAGAAG
 GTCTGGTCCATGCCTCAGCAGACCAGACCCGAACCCACACCTGATGGCAGACACAGATATCACTGAGAGA
 AAGAGCAGTCATCAAACGCCAAGTATGGTCCAGAAATGAATTGATCGAGTAAATTAAGAGGTGGAAGAA
 TTGGCAGAAATTTCAAATACACGAATGTCTGTTTTAGAAGAAACCACAAGTGCCTCAAGTTAGGTACAG
 ACAAGAAGCATCCTTTCCCTTCCTCCGTATCATCTGCTCTTCTGGTGGTCTGCTTACTGGGAAGACTGCT
 CCTCCTCTCCCTAGAGCGTTACCTGTGAGTCAGCCTTGGCAAAGAGCAGCAAAAGGAAGCTGAGGGGC
 AGGAATGTGTTCCCTTCTCCTGTGGGGCCTCAGCCAAGAAAACAGGAAGTGGTGTGTTGTATTGCCCT
 GCGCTTTCTTACCTATTACGTCTGTACAGAGCCCCAAGAGAAAAATGCAATTTGGTAGTGCACTGATAAC
 TACATTATGAATCATTGGTCCAAGAATGATAAAATCCAAACAATGCTAGGGGAAGGAGATGGAAGAG
 AGGCAGAACAAAACCAACACAATCTCATTGTGAAGGACTGACAATTACAATTAAGTTAATGTATAACT
 TCCAAACGGAAAATTTT

OTI All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence [imiRBase](#). The
Disclaimer: sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should
 affect the function of the mature miRNA.

E. coli Kanamycin (25 ug/mL)

Selection:

Mammalian Neomycin

Cell

Selection:

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Both mature and minor (star) miRNA can be over expressed from the construct.

Contains:

RefSeq: [MI0021954](#)



[View online »](#)