

Product datasheet for **SC402767**

MIR669a-6 Mouse MicroRNA Expression Plasmid (MI0014058)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	UGUGCAUGUGUGUAUAGUUGUGUGCAUGUUAUUAUUGAAUUAACAUAACAUAACACACACAGUAUAAACGCAAGCACA
Mature Sequence:	AGUUGUGUGUGCAUGUUAUGUCU
Full insert:	>MI0014058 AGTGAAGTTCGTGGTAAAGGCATGATGAGTTAGAGGTCATCCTGGCATTTCATGTGATTGCCCTATGTA TGTGTGTATGTGTTGATGTGTGTGTACATGTACATATGTGAATATGATATACATATACATACACGCACAC ATAAGACACATATGAGCACACAGACACATGGGAATGGCACTCATCAATCCATTGAGTGATAAGTGACTT AGAGTGAAATGCTTGCTAGAGTTAAGAGTATCAGGAGTTCAAATTCATACTTATCCTCATGAGAGATCA GCTTGTGATTCTCCATGTATGTGCATGTGTGTATAGTTGTGTGTGCATGTTTCATGTCTATATTTGAATA TACATAACATACACACACGTATAAACGCAAGCACACATACACAGACACAGGAATGGCATTGCGCAATC CCTTTGAATGATAAGTGATGTAGAGTTTATAGAGATTATGCAATAGGGGGCATTGAATGGCCTCACTCACC AGGTGTAAAAGATTGACAAAGTCAGGAGGATAAGGATTGCAAATTCATGCTTAACTTCATTTGGAGACAT CTTGGTATTAATGTCATGCAACCCATTTATGTGCATGTGTGCTTACTTGTATGTGCATTGTATGTGTTGA CATGATATACATAACATAAACACACACGTGTATCTGCAAGCACACACACAGGAAAAA
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0014058
Synonyms:	mmu-mir-669a-6
Locus ID:	100526485


[View online »](#)