

Product datasheet for **SC402783**

MIR6903 Mouse MicroRNA Expression Plasmid (MI0022750)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMVIR)
Precursor:	CCCUUUUGGAAAAGCAGUACUUACCACAUUGUCUUCAGGUGGCACGCGUGUGUUCACUGUCGUGUGUCUGUUUUUCCAGGUAUCA
Mature Sequence:	UUGGAAAAGCAGUACUUACCA
Full insert:	<pre>>MI0022750 ACAGCTATCTAGAGCTGACAGTAGGAGTGTTAGTCCCTAGCTGGTAAGGCTGAAGTGGTTTCCTTGAAC CAGATGTCTGAAGTGGTATCCTTGAACCTCAGATGTCAAAGCCATCTCATGCTTTCCAGCTACCTAGCAA CTGTGTAAGGCACAAGACCAAGGCTACCTTTCTCGGTTTCAGAGTGAACCTCAGCGAATCAATTGCCTGTGG GTGGTTACTGTGCGCCAGCAGCAAAAAGACAAAGGAAGGTGAGAGAGTAGAGGCAGAATAGTAGAGACA GGGAGAAGACCAAGGTGTCCTTTTGGGAAAAGCAGTACTTACCACATGTCTTCAGGTGGCAGCGGTG TGTTCACTGTCGTGTGTGTCTGTTTTTCCAGGTATCATTGACAGCACTGCTGGGGAACAGCGGCAGGTA GTGGCTGTCAACCGGTGATGGAACAAATGATGGACCTGCTCTGAAGAAAGCAGATGTCGGGTTTGCTATGG TAAGCTCCCCTTGGGTCACCTTCTCCATCAGGAGCCAAGTGTCTGCTTCCGATGCCTTCGAATCCTGCCT GTGCTCCTCCACCTGCACCCACCCCCACCTCTGTCTCTACCGACCTGGCCAGTAGTGTGTGAAAA CCAGACGATGGTTGTTTCAGGGAATCACAGAGTGGCTTACCTGGAAAAACAAGAGCACT</pre>
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0022750


[View online »](#)