

## Product datasheet for **SC402818**

### MIR6937 Mouse MicroRNA Expression Plasmid (MI0022784)

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | MicroRNA Expression Plasmids                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Species:                  | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Vector:                   | pCMV-MIR (PCMVMIR)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Precursor:                | UAGCUGUAAGGGCUGGGUCUGUGUAUCGGUCUCCAGGAAACUCACAGGUUCCCCUUGCACUAG                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Mature Sequence:          | UAGCUGUAAGGGCUGGGUCUGUGU                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Full insert:              | <p>&gt;MI0022784 TCACTGGGACCATTGAGCTTAGTGGCATAGTTATTTAGAGTGGAGCATGGTCTCAGGTAGGTAATCTCA<br/> GCTTCTGGTAGGTGAAGACAGGAGGAGCAGAAGTTCAAGGTCATCTCAGCTACTTAGCAAGTTTGAGGAC<br/> AGCCTGTGCTACATGAGACCTTGCCCTCAGTTTGCCCTACCACTTGGGTTCTGTGGATTAATTTAAGAC<br/> ATCATCTCGGTGGCAGGTGCCTTTACCGACAGTGTCATTTCTGGCTTGGCTCTTGACACTTTTAATACC<br/> AAGAATTATTCAGTGTAGTCTAGCTGTAAGGGCTGGGTCTGTGTGATCGGTCTCCAGGAACTCACAGGT<br/> TCCCCCTTGACTAGATCAAGATGTTCTTTGAGGAACATCTGCATCTGGATGAGGAGATCCGCTACATCC<br/> TGGAGGGTAGTGGTACTTCGATGTGAGGACAAGGAGGACAAGTGGATCCGGATTTCCATGGAGAAGGG<br/> GGACATGATTACTCTTCCTGCCGGCATCTATCACCGCTTCACACTGGACGAGAAGGTGAGAGTGTACCAT<br/> GGTGTGGTCTGTGCATAGACATCGCAGGTGACCATGATGACCCACTAGGTGTGTTGTGTGAGGCTCATCG<br/> CCTAGATACATGGCCATTACAGCTCTTCAGAGA</p> |
| OTI Disclaimer:           | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Components:               | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Also Contains:            | Both mature and minor (star) miRNA can be over expressed from the construct.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| RefSeq:                   | <a href="#">MI0022784</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |



[View online »](#)