

Product datasheet for **SC402819**

MIR6938 Mouse MicroRNA Expression Plasmid (MI0022785)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Mouse

Vector: pCMV-MIR (PCMVIR)

Precursor: AGUGUUUGGAGCUAGCUCCAGAGAGAGUGAAGUCAUGCUUGUCAUCUGGGGCGUGUCUCCUAG

Mature Sequence: UUGGAGCUAGCUCCAGAGAGAGU

Full insert: >MI0022785 TGCTTCTCAGAGCGGCATGCAATTTAAGTTTTATGGCCAGGGACATAGCTCAGCGGTAGACCATATGCTT
AGCATGTGCAAGGCCCTGCTCAGTCCCTAAAATGAAGTTTACTTCTGGAATTTCCATTATTTTTTTG
GATCAGAGTCAACCACTAGTAAGTCTTCTAGAAAGTGAAGTGTGAATGAGAATTACTGTGATTGGA
AAGGAAATTTGATAATAAAATCCTAACACATATAAGCATCAGGAAAACAATACGACTTCCAAAGGTCATA
GCAGCTTTTGAAATGATGGAGTGTGGAGCTAGCTCCAGAGAGAGTGAAGTCATGCTTGTCATCTGGG
GCTGTCTCCTTAGGTCTTACCACAAACATTTAGAAGATCTGGGGCTTGAGTTGCTCTTCCAGATGCTAG
TGATTCTCTGATCCTGGTGGGAAAAGTGCCGCTCTGCTTTGTAGAGAGAGAAGCTAGTGAGCTTCGAAGA
GGACGCTCTACTGTGACTAAGAGTATTGTGGAGGTAAGACACGTCCTCCGCGTTTGCACTCTTGCTGAGG
CCGCTAACTGCAATCTTTGTGTCGGTCTGCACATGTTTGTGTCATGTTACAGTGCCAGATGATGGAAA
TTGACAGATGTATGAGGTCACAGTAAGGAGGAA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0022785](#)



[View online »](#)