

Product datasheet for **SC402825**

MIR6944 Mouse MicroRNA Expression Plasmid (MI0022791)

Product data:

- Product Type:** MicroRNA Expression Plasmids
- Species:** Mouse
- Vector:** pCMV-MIR (PCMVMIR)
- Precursor:** GUGAGAGCGGGGGAGUGGCAAGUGAGACGGAUGGGCGGGGCGGUCAGACUGCUAGCCCCGCCAGCUCUCUAAACUCUCCCUUGUGCCUCAG
- Mature Sequence:** GUGAGAGCGGGGGAGUGGCAAG
- Full insert:** >MI0022791 TGGTGTGACTTAAGCTAGCTAAGGGACCAGGAATATTCTGAGTGATCCCAGATACCTGGTTAGATGGGACGTAGATGAATCAGCGGCTCAAGGCCTGGCTGGCCGGTCTGGGGATGAGGTGGCTGCAGACATTGAAAGAGAGGCTGTGGACCTGAATTCTGAGAGGAACAAAGCGCCACCTTTGAGATTCTGCGGGGTTGAAAGAACCCTCAGGACATTGTCCCTTTGGTTTTCCAGGGTGACCAGCGGGATGAAGAGGAAGAGTCCCAGATGAAGAAGTCCGAGTCAGAGGTGAGAGTGGAGTGGCAAGTGGGCGGGGCTGGCTCAGACTGCTAGCCCCGCCAGCTCTCTAACTTCCCTTGTGCCCTCAGGAGCGGCTGCCATCATTGCCCAGCGGCCTGATAACCCACGGGAGTTCTTCAGACAGCAGGAACGAGTGGCATCGGCCTCTGGTGGCAGCTGTGACGCGCCTGCGCCTGCACCTTCAACCACCGACCAGGTAGTCTTGGCCCCCTCCTTGCCCAATCAGTACTCCTCAGACTATGGCCCTAATCCCTCTGGGACTTCTTCTGGAAAAGAAGCCTTGCCCAAGGAAGGACTCCTTGGAGAAGGCAAGACCCCTTCAGTGTGCTCCACCCAGCTTCCAGGTGGTGGACCC
- OTI Disclaimer:** All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
- E. coli Selection:** Kanamycin (25 ug/mL)
- Mammalian Cell Selection:** Neomycin
- Components:** 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
- Also Contains:** Both mature and minor (star) miRNA can be over expressed from the construct.
- RefSeq:** [MI0022791](#)



[View online »](#)