

## Product datasheet for **SC402837**

### MIR6956 Mouse MicroRNA Expression Plasmid (MI0022803)

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | MicroRNA Expression Plasmids                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Species:                  | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Vector:                   | pCMV-MIR (PCMVIR)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Precursor:                | GCAGGCAGUGAGUGGGUGGGACUGGCCUGGAGUCAGUCACCCUGACCGGCCUAUCCUCUCAG                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Mature Sequence:          | UGAGUGGGUGGGACUGCCU                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Full insert:              | <p>&gt;MI0022803 CAGCCACTGTAGCCTGGGTTCAAACCACCATGTATCTAAACCCAGTAATGTTGCGATGGGACTCTGGGC</p> <p>CAAGCCGTGGCATGAGCGCCCATGGTTATCTGCCCCTGTCCCGGCTGGCTAGGCTGGACATGGAGGC</p> <p>GCTAAAGCAGGAGCTGCAGGTGTATCTGAGAGGTAAGTGCCTGGAGATCGGGGCCCTGACA</p> <p>AGGCAGGCAGAAGAGCGTGAGCACACCTGCGGCGCTGCCAGCAAGAGGGCCAGGAGCTGCTGCGCCACA</p> <p>ACCAGGTAGGCGGGGCCCTGGCAGGCAGTGAGTGGTGGGACTGGCCTGGAGTCAGTCACCCTGACCGGC</p> <p>CTATCCTCTCAGGAGTGCACTCCACCTGTCCGAGGAGATCGACCGTCTGCGCAGCTTCATCGCCTCGC</p> <p>AGGGCACAGGCAACAGCTGCGGGCGCAGCAATGAGAGGAGTTCCTGTGAGCTGGAGGTAACGTCCTCCC</p> <p>TCCCCCTCCTCCCTCATGCCAGGCGACCCCTGGCTTGAAGACCTGTGAGGTAGCTCTCTAGAAGGCTGGC</p> <p>TGTCACCAAAAGGTGTCACTCAAACCTGAGGTCTTCCAGGCTAAGCCTTCTAGCCTTTGATCCCGCTC</p> <p>TGTCATCTGCTTGGGAGTGCTGCTCAGTC</p> |
| OTI Disclaimer:           | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Components:               | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Also Contains:            | Both mature and minor (star) miRNA can be over expressed from the construct.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| RefSeq:                   | <a href="#">MI0022803</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |



[View online »](#)