

## Product datasheet for **SC402860**

### MIR6978 Mouse MicroRNA Expression Plasmid (MI0022826)

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | MicroRNA Expression Plasmids                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Species:                  | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Vector:                   | pCMV-MIR (PCMVIR)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Precursor:                | AAAGGAAGGGGUGAGAGAGAAGCUGGUGCCAUCCACGGCUUCACUCACCCUGCAG                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Mature Sequence:          | AAGGGGUGAGAGAGAAGCUGGUG                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Full insert:              | <p>&gt;MI0022826 TGCACACATACACATGAGTATAAGGTAGGCATTTGTATCACCTTCTCCTTAGAAGAAGTTGAGACATAGA<br/> GAGGTTAAGTAACTGGTCTCCAGTCATACAGGCAGAGAGGCAGGGCCTTTGGAAGGACTTTCTCTGAAG<br/> AGTACCTGCTGACCATCTGTGTATCTGTTGTCAATTTATCCACCACCTCCACAGCCAGAAGGAATGAG<br/> AATCAACAAAAGTGTGGCCATCCATACACTGGACCCAGAGAAGCTCGGTCAAGGTGAGGCAGGCAGAGAG<br/> GGCAAGGTGGGTAGGGTACAAAGGAAGGGTGAGAGAGAAGCTGGTGCCATCCACGGCTTCACTCTCAC<br/> CCTGCAGGGGAGTGCAGAAGGTGGATGTGCCTGCCGACAGCCTAGCGACCAAGTGCCAGACACAGACT<br/> CTGAGACCAGAATTATCCTGCAAGGTAAGATGGCCTCTCCCTGCCATGCAGGGGTGAGGGTAGCTGGGA<br/> GCCAGAAGTTAACTGAGCTCCATCTCATGGGCTGAGGGTGAGGAGAGAAGAGAATAAGGAAAGGGAGCGG<br/> GGGAAGACCCTGGTCCAGGATGCCAGATTCAAATCAATTTCAAGTTGATTAACCTGGGACTCATTGTGTT<br/> AACATACTGGCACCCCTGTGCCTCAGTT</p> |
| OTI Disclaimer:           | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in <a href="#">miRBase</a> . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Components:               | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Also Contains:            | Both mature and minor (star) miRNA can be over expressed from the construct.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| RefSeq:                   | <a href="#">MI0022826</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |



[View online »](#)