

Product datasheet for **SC402914**

MIR7032 Mouse MicroRNA Expression Plasmid (MI0022881)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Mouse
Vector:	pCMV-MIR (PCMV MIR)
Precursor:	CGUGGCCAGGGUGGUCCCCGAGAGGGUGGUCCCCACAUGGGCAAACCAGGGGUCCAUCCUCUCGGUACCGCCUCGAG
Mature Sequence:	CCCAGGGUGGUCCCCGAGAGGGUGGU
Full insert:	>MI0022881 ACCCCCCAATGGCTCCATGAAGGATGACAGGTGAGTCCCGAGGAGCCAGCTGGGGTGGGGAGGCTAGGAA AGGCCTTCTGTCTGGGAAGGTCATAGGTAGGACCTCCAGAGCAATGACACACCCAGCCCTGAACCACTGC TAGTGACATCTGAGATCCGGTGTGGGCAGTCAGGAACCTTCTGGATGATGAACTTGAGCCAGGATCTC TGTCCCAGGGACCTCCAAATCGAGAATTGAAGAGAGAGGTGGAGACCTCCGTGCTGAGCTGGAGAAG ATTAAGATGGAGGTGAGGGCGTGGCCAGGGTGGTCCCCGAGAGGGTGGTCCCCACATGGGCAAACCAG GGTCCATCCTCTCGGTACCGCCCTGCAGGCACAGCGGTACATCTCCAGCTGAAGGGCCAGGTGAATGG CCTGGAGGCAGAGCTGGAGGAGCAGCGCAAGCAGAAGCAGAAGGCCCTGGTGGACAACGAGCAGCTGCGC CAGGAGCTGGCCAGCTCAAGGCCCTGCAGCTGGAGGGCGCCCGCAACCAGGGCCTTCGAGAGGAAGCAG AGAGTGTGTACCTCTGGGGCGGCACCCAGCAACAGGGTGAGGGTGGCTTACCCCTGCCTGCCTACTCC ACTCAGCCTCACCTCTCCGACAGGAAGGCCAGTGCCACGGAGGCACGCT
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0022881



[View online »](#)