

Product datasheet for **SC402915**

MIR7033 Mouse MicroRNA Expression Plasmid (MI0022882)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVIR)
Precursor: UUUAGUCUCCAGGAGUCUGAGGGGCAGGGUCGUCUCCAGCAUCUUGGGCUUGGCCUAACCUGCCUCAUUCUCUUGGGCAG
Mature Sequence: UCUCCAGGAGUCUGAGGGGCAGG

Full insert:
>MI0022882 AGTCCCAGGAGAAAAGGGTAGGTAGAGAAGCCATGGGCCCGAGAGGAGCATCTCCAGTGAAGTTCACCG
TGTTACTCTGGGTGGGCAGGAAGATCGTCTTGCTCTGAGTAGGGAGGTGAGCAAGCCTGTGTTTCTG
TCTGTAGTTACAGATAACAACCTGGGAAGCAGCAGAGCTCATGGACTTCCCCACCCACCTGGCAGCT
CCGGGCAACGCAAACGCAAGATCCAGCTGCTGCCCTCCCGGCGAGGGGACCAGCTCACCTTGGTATGTCT
TGTCTATCCTTCTCGGCTCTTTAGTCTCCAGGAGTCTGAGGGGCAGGGTCTGCTCCAGCATCTTGGGCT
TGGCCTAACCTGCCTCATTCTCTTGGGCAGCCTCCACCCCAGAGCTTGGCTATTCCATCACTGCTGAA
GACCTGGACATGGAGAGGAAAGCCTCACTGCAATGGTTAACAAGGTCTTGGAGGACAAGCCAGGTCAGG
CTGGGCGCAGTGTGCTCCGAGTGTCTGTATTCTTTAGCTTTCTTGAAGAGCAGTTGTAGAATA
CTTGATTTTCCATTCATTATCCAAGGCTGTCTCAGGGTCATAATCCTCCTCCTGCTTCAGCCTCCAAA
GTGCAGGGTTGCAATGTGCACCATTACAACCAGCTCTTAGATGTTTGTCC

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0022882](#)



[View online »](#)