

Product datasheet for **SC402958**

MIR7075 Mouse MicroRNA Expression Plasmid (MI0022925)

Product data:

| | |
|---------------------------|--|
| Product Type: | MicroRNA Expression Plasmids |
| Species: | Mouse |
| Vector: | pCMV-MIR (PCMVMIR) |
| Precursor: | GGCCGCGGGGAGGAGGACAUGUUUUUGGGGGCAUGACUUGGGGUCCGGGGCUCGCGGAGCCAACCAUGUCUUCUUUCCAG |
| Mature Sequence: | CGGGGAGGAGGACAUGUUUUU |
| Full insert: | >MI0022925 AGGGCGGGGAGGGCGAGCGGCGGTGCAACACGATGGCGGGAGGAGGGCGGAGCAGCGATAGCAGCGGGCG GGCGGGGAGCCGACGGGCATCGCGCAGCGGGCGGGCGGCTCGACGGGGGCGACACGAGCCAGGGTTGGGA GGCGGGCCGAGCGGGGCGCGGGGAAGCTCGCCTGGAGGAGGCGGTCAACCGCTGGGTGCTCAAGTTCT ATTTCCACGAGGCGCTGCGGGCCTTTCCGAGTAGCCGGTACCGGGACTTCAGGCAGATCCGGGACATCAT GCAGGGTGAGAGCCGGCTCGGGCCGCGGGGAGGAGGACATGGTTTTGGGGGCATGACTTGGGGTCCGGG GCTCGCGGAGCCAACCATGTCTTCTTTCCAGCGTTGCTTGTGTCAGGCCCTTGGGAAGGAGCATAACGGT TCCCGGTTGCTGCGGGTTATGCAGTGTCTGTGCGCATTGAAGAAGGAGAAAAATTTAGGTACGGCTATTT TTTGTGTGTTCTTTCTGGGCTTACTTACGTAGGGTTTTGTTATTTTCAAACCTAACATCCGGAAGTAA AGATTCCGGGTTAAAGTGCACACTACTGGTTCACCGCGGAGGGCTTTTTTCATTCAACAGTGGTTGAATTCT CTAATCCGGTTTTCTCCAGTTAAAGTTTCTTGCTGTCTTTTGTGTTTTGTTT |
| OTI Disclaimer: | All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA. |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Mammalian Cell Selection: | Neomycin |
| Components: | 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube |
| Also Contains: | Both mature and minor (star) miRNA can be over expressed from the construct. |
| RefSeq: | MI0022925 |



[View online »](#)