

Product datasheet for **SC403056**

MIR7681 Mouse MicroRNA Expression Plasmid (MI0025029)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Mouse
Vector: pCMV-MIR (PCMVIR)
Precursor: AUCCUGUCCUUGCCCUCUCUGUUGGCAUCAGAAACCUACACCUACACUACAGAAAGGGCACUGACAGGAUAG
Mature Sequence: AUCCUGUCCUUGCCCUCUCU

Full insert: >MI0025029 TCCTTTCTCCCATTTATTTTCATCATTTGATCATGTATTTCAGTCACTTAAAATTACTGTTTCATGCCGGGCG
 TGGTGGTGCAGGCCTTTAATCCCAGCACTTGGGAGGCAGAGGCAGGTGGATTTCTGAGTTCGAGGACAGC
 CTGGTCTACAAAGTGAGTTCAAGGACAGCCAGGGCTATACAGAGAAACCTGTCTTGAAAAACAAAAAC
 CAAAAAAAAAAAAAAAAAAATTACTGTTCAATATTTCTTCATTTTGTAGCTGAAATTATTTCCACTTTG
 GTTACTAGTTTCTGGGCTCTATCCTGTCCTTGCCTCTCTGTTGGCATCAGAAACCTACACCTACACTAC
 AGAAAGGGCACTGACAGGATAGAGCCTAGAAACCCAGTGTCAACAACCTGTGGGTAGAGGGAGTAGACAA
 ATGGATTGTTTCTGATAATTGAAACCTTCCCAGAACATCCTCTCTCCTTTTTTCTTTTGGTGTACAG
 GAAATAGACAAGTGACAGGGGTAGGAAGCACTTGATCACTAGTGGTCATCAACAAGAAGGAATGTAGCT
 CTGCTCAGAAGTGATGCTATCAGATCAGAACATGGCAGCCCACTTCAGGATGATTGCACTCAACAACGT
 CTGTGGCTACTCAAATGTATCCACTGCTAGGGCTTCTGAAGT

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](https://www.ncbi.nlm.nih.gov/mi/). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0025029](https://www.ncbi.nlm.nih.gov/nuccore/MI0025029)



[View online »](#)