

Product datasheet for **SC403108**

MIR100 Rat MicroRNA Expression Plasmid (MI0000885)

Product data:

Product Type: MicroRNA Expression Plasmids
Species: Rat
Vector: pCMV-MIR (PCMVMIR)
Precursor: CCUGUUGCCACAAACCCGUGAGAUCCGAACUUGUGCUGACCAUGCACACAAGCUUGUGUCUAUAGGUAUGUGUCUGUUAGG
Mature Sequence: AACCCGUGAGAUCCGAACUUGUG
Full insert: >MI0000885 TCAGAATCGGTCTCGTTCTGCCTATAGGCCATACTGACATTAGTTTAGCCTACATAGATGAGGATTCCGA
TAAATGAAGCTCACCTGAAGCTCACTTCGTCTGTGAGTGAAAAGTTAGGTCTACCAGAGGGAACAGGAG
ACTTCTCTTTAAAACCTTTCTGTCTTTGGCATGAAATGGCCCAATGGCAAAATCCGTTGTTGGACTAAA
TGAGCCAAGGACAACTTCATTGTTGGACTAATGGTTTTAACGTTTTATCCTCAGACAAGTCAACCCCA
CACCAGAGAAGATACTGAAGCCTGTTGCCACAAACCCGTAGATCCGAACCTGTGCTGACCATGCACAAA
GCTTGTGTCTATAGGTATGTGTCTGTTAGGCGATCCCACAGTCCAGGGGCTTTGCCTGTTGCCTCTAGT
ATCTGCTCCTCACCTCCTGCCCTCCCTGTCTCCCTGTGGTTGCCTTTTCAGGGTTGTAGAAACGGAAAGAG
GGCTGAAGAGGAGGGAACGGAGCTCTTGGGCCACGATGGATAAGGTTATTTAAGAGACGCACTCTCATA
TGCTTCTAACTGGTCTCAGCTGTCAGGTTAGGGGGCTGGGCTGCCATGAGGCAGCAAGAATTCAGGGA
TTCAGACAGCTGACTGCTGACTGTGGAGTAAATAGGCTCACACAAAACAA

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)
Mammalian Cell Selection: Neomycin
Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq: [MI0000885](#)
Synonyms: rno-mir-100
Locus ID: 100314234



[View online »](#)