

Product datasheet for SC403149

MIR137 Rat MicroRNA Expression Plasmid (MI0000910)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Rat
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GGCCUCUGACUCUCUUCGGUGACGGGUUUUCUUGGGUGGAUAAUACGGAUUACGUUGUUUUGCUUAAGAAUACGCGUAGUCGAGGAGAGUACCAGCGGCA
Mature Sequence:	ACGGGUUUUCUUGGGUGGAUAA
Full insert:	>MI0000910 TGGCACCTAGGCTCTGCTTGGCGGAAAGAAAGCACTTGCTGTTGTGTGTTAAGTGGCTGTTCTGTGAT TTGTCTTTCAGAATTGGAAATAGAACAGCCATTTGGAAGATCTGGGCAGGAAGCAACAGAGCACAGCTTT GGAGCCTTCTTAGGGAAATCGAGTTATGGATTTATGGTCCCGGTCAAGCTCAGCCAGCCCTAGCCAGG GGCGGGCTCCGACGCAAGAGTCTTTCTGGTGGTGGCGGCGGCAGCGGCAGCAGTGACAGCGGTAGCAGA GGCAGAGGCAGCGGCAGCTTGGCCCTCTGACTCTCTTCGGTGACGGGTATTCTTGGGTGATAATACGGA TTACGTTGTTATTGCTTAAGAATACGCGTAGTCGAGGAGAGTACCAGCGGCAGCGGGCAGCAGCCGCC CTCCCAGCCCACCAGCTGGCCACTAAAAGCTCCCGTTGCCAAGGTAGCACTTTTTATGTTTTCTTTTCC TCGGGTGGTTTTCTGTTCTGGATCTTCTGGCACCTCTGGTGACCCAGACCCTTTTAGAAGGAAAGGGGA GGGTGCGCGTCCGAGTTGGGGCGGGGGCAAGCGGAGGAGAAGGGAGGGGAGAGCGGAGAAAGGAA GAGGTCTCTGAAGCCCGGGCTGTGGGGCCATTTGAGTAGTCTCTCTCTCTCTCTCTCTCTCTCTCTCT TC
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000910
Synonyms:	rno-mir-137



[View online »](#)

Locus ID: 100314031