

Product datasheet for **SC403153**

MIR140 Rat MicroRNA Expression Plasmid (MI0000611)

Product data:

Product	MicroRNA Expression Plasmids
Type:	
Species:	Rat
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GUGUCUCUCUCUGUGUCCAGUGGUUUUACCCUAUGGUAGGUUACAUCAUGCUGUUCUACCACAGGGUAGAACCACGGACAGGAUACUGGAGACC
Mature Sequence:	CAGUGGUUUUACCCUAUGGUAG
Full insert:	>MI0000611 AGCAATGTGCTGCCCTCACCTATACAAGTTAACAGCTCTGTGGAAAGCATCACCTCGTACAACAGTCA GAGGGTTCTAGCAGTGGGAAGGAGCTTGTGGCTGTGATAACCTCTTATCCTGCTGACCTTGTGGTGCGG GGATATCTTCTGGCTTGGTGTGGGCTCCCCGCCTCTGATCCTATTTGCGGTGGTGTAGTCTTCTGTGT TGATCCCATCCTGCAGATGAATTGTTCTTTTTCCGTGGTGACCTCTCCAGGCTGTGCTTGGTGGGCATC TGGTGTGGCTCCCGCCTGTGTCTCTCTGTGTCTGCTGCCAGTGGTTTTACCCTATGGTAGGTTACAT CATGCTGTTCTACCACAGGGTAGAACCACGGACAGGATACTGGAGCACCTCTGCATCGAAGGACTCCAT CTACACAGCCACACACAGCCCAGCAGCCAGGGGCCGGCTCCCTTGACCTCTGTTCCCTGACTGTTTAGCT CAAATCTGAAGGACTCTGGATTCTGGCTTGGATTCTATGTAGGATGCAGATAGATGCTAAGGCGTCT GTTGGGTGGGTGACAGTTGAGATGAATGGATCCAAGTTTGGCTTAGGGTACGCTAATAGGGTCTAAGGTG AAGGGAGGGCAGCGTCAAGTGAGGGACCTGTGGGTAGCGTGAAGGAATGAGGTTTCTAGAGGCCCAA
OTI	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in
Disclaimer:	miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000611
Synonyms:	rno-mir-140
Locus ID:	100314276



[View online »](#)