

Product datasheet for **SC403203**

MIR19b-2 Rat MicroRNA Expression Plasmid (MI0000848)

Product data:

Product Type: MicroRNA Expression Plasmids

Species: Rat

Vector: pCMV-MIR (PCMVMIR)

Precursor: ACAUUGCUACUJACGGUUGUUUGCAGAUUUGCAGUUCAGCGUAUAUGUGGAUUAUUGGCUGUGCAAUCCAUGCAAACUGAUUGUGAUGAUGU

Mature Sequence: AGUUUUGCAGAUUUGCAGUUCAG

Full insert: >MI0000848 GAAGCAGCTTATAATCTACTGCCCTAAATGCCCTTCTCGCACAGGCTACTGCTTACTCCGCATTTGAAA
AGCATATCCTGAGAAATGATTTGAATGGGACAAGTCTTAAAGGCTTTTTCACTGTTTCACTGAAAAACA
AGATGGTGGGGACTGTGCAGTGTACAAGGACAAGCTCCGATCCTAGTAGTGCCAAAGTGCTCATAGTGCA
GGTAGGTTTTGCTGCACTCTACTGCAGTGTGAGCACTTCTGGTACTCTTGAATCAACCAACCCTCGTT
GACGCAACTACTGTGCGAGGACATTGCTACTTACGGTTAGTTTTGCAGATTTGCAGTTCAGCGTATATGT
GGATATATGGCTGTGCAAATCCATGCAAACTGATTGTGATGATGTGTGCCTTTCTGCATCTTTATGGAC
ACCTTTAGCACCATGCCCATTCACACAGGTGGGATTAGTGCCATTACTTGTGTTAGATAAAAAGTA
TTGCACTTGTCCCGCCTGAGGAAGAAAAGAGGGTTTTTAATCATCTTTTAGTTTCTGAGAGTTGTAAGA
TATGTTAGACTATAATTATGAAATGAACTGTTTTGCTGTTATCGGGTGGATCACGATGCAATTTTGATTA
GAGTAATGGGAGAGAAATTGCACGGTATCCATCTGTAACC CGCAGGACTTTTGGTGACTGCAGTCG

OTI Disclaimer: All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in [miRBase](#). The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.

E. coli Selection: Kanamycin (25 ug/mL)

Mammalian Cell Selection: Neomycin

Components: 1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube

Also Contains: Both mature and minor (star) miRNA can be over expressed from the construct.

RefSeq: [MI0000848](#)

Synonyms: Mir19b-2 rno-mir-19b-2

Locus ID: 100313999



[View online »](#)