

Product datasheet for **SC403210**

MIR203b Rat MicroRNA Expression Plasmid (MI0015399)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Rat
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GCUGUCGCCGUCGCCGUCGCCGGGUCUAGUGGUCCUAAACAUUUCACAAUUGCGCUACAGAACUGUUGAACUGUUAAGAACCACUGGACCAGGCGGCCUGC
Mature Sequence:	AGUGGUCCUAAACAUUUCAC
Full insert:	>MI0015399 CCCTGGAGCAAAGAACACCCCTTGACCCAGAGGCTGCCGCGCTTCCTGCCAAAGATTCAGTGGTTTTG GACTAGGGAGGATCGAGCAGGAGAGAAGCAAGGTGACGCTAGCGCTTGGCCAGCGCACGTGTCAGCAG TGTGGCGAAAGAGCAGGCTGCCGAGGTCCCGAGCCAGCGTGCCTCGGATCGTCTGCACTCCA GCCAGGCAGCGAGAACGCACGGGTCTCAGCGATCCGGGTGCCTTCTCCCTTTCAAGCCGACCTTGCC GTCGGGGCCGCCCGTCTGTTGCTGTCGCCGCTGCCGTGCGCGCCGGGTCTAGTGGTCTAAACATTTAC AATTGCGCTACAGAAGTGTGAAGTGAAGTAAAGAACACTGGACCAGGCGCGCTGCACACAGAAGCTCTGC CCAGCCCGGAACCCCGGGCCACGGCTCAGCGTCTCCTGCGACCCGGCCCGCTGGCCCGGCTGCACTC CGTGGCGCGCGCCGGCGTGCCTGCGCTTTTATGTGGCACGGCCGCGCGCTCGGCGCCGGAACGCGC CAGCCGCGCTGGCCGGACCGGTTTTTCCAGGCCCGCGGGTCCGCGCGGGTTCGGCCCGAGGCCCCAGGA GCTGCTCCGCGTCCCGCGCGGGGTGGCCCTGAGTACCCGGCCGACCCGGCGCCAGCTCGACGCATC TGGCCGGGTTCGGCG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This may affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0015399
Synonyms:	rno-mir-3545



[View online »](#)

Locus ID: 100526564