

Product datasheet for **SC403219**

MIR21 Rat MicroRNA Expression Plasmid (MI0000850)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Rat
Vector:	pCMV-MIR (PCMVIR)
Precursor:	UGUACCACCUUGUCGGGUAGCUUAUCAGACUGAUGUUGACUGUUGAAUCUCAUGGCAACAGCAGUCGAUGGGCUGUCUGACAUUUUGGUAUC
Mature Sequence:	UAGCUUAUCAGACUGAUGUUGA
Full insert:	<pre>>MI0000850 TATAAATCAAAGGACACTAGCAAGAAAGGCTCTCGTGTGTGTTAGTTTTGCTGTTTGTACTCCGGCTTCA TCAGGCGGAAGTTTACTTCTCTGTGTTTCTGTATTTTATTCATTGTGATTTTCCTTTGACATTTTTCTG ATTGACCTGTTCATTTTGTGTTTGTGTTGGTAGGAAATAAACCGCTTTCTTTCTAGAAATTGGCATTAAAG CCCCGGCAAGCCAGAGACGTTTGTGTTTAAACCCCGCTGAGCATCTCGTGCTTGCCCCCTGTCT GGTCGTGACATCTGCATGGCTGTACCACCTGTGCGGTAGCTTATCAGACTGATGTTGACTGTTGAATCT CATGGCAACAGCAGTCGATGGGCTGTCTGACATTTGGTATCTTTCATCTGACCATCCATATCAAATGTT TTCATTCAAACATTACCCAGCATCATGGTTTGCAGTCAGAAATTGGTCCTTATGTCTCATGAGACTTG CAGCCCTTATACTGTGTTACATCAGCAACATCGAGAGGAACTTCCTAGGAAGTTGGGGTAGTCGTCAGA GTCCTAACTGAATAAATAAAATGAAAAGTCACCTCGCCTTGTGAATTTGGAAAGTTGTCCTTCGTTGG TTTTGGTTGGCTGGTTGGTTTATGTTTTTTCTTTACAGGGACTAACATTTGTTACCTGAC</pre>
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence in miRBase . The flanking sequence of a miRNA clone may differ from the NCBI reference with respect to biological polymorphisms. This should not affect the function of the mature miRNA.
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0000850
Synonyms:	rno-mir-21
Locus ID:	100314000



[View online »](#)