

Product datasheet for **SC403224**

MIR215 Rat MicroRNA Expression Plasmid (MI0003482)

Product data:

Product Type:	MicroRNA Expression Plasmids
Species:	Rat
Vector:	pCMV-MIR (PCMVMIR)
Precursor:	GGUGUACAGGACAAUGACCUAUGAUUUUGACAGACAGUGGGCUGCGUGUGUCUGUCAUUCUGUAGGCCAAUUAUUCUGUAUGUCUCUCCUCCUACAA
Mature Sequence:	AUGACCUAUGAUUUUGACAGACA
Full insert:	>MI0003482 AGAAGCATACACACTGAGTGACTAGTGCTTTAAAGCCAAAGCAGTTCTCCGTCAAGCTTATGGAGTCATC ACGTGTAACAGCAACTCCATGTGGACTGTGCACAGATCCCAGTGGAGCTGCTGTACTTTTGGATGGCCTC CAGCCTGCTACAAACTGGAGAGGAAGANNNAGAGTGTAAATGTTTTAAATCCAAG
OTI Disclaimer:	All miRNA clones were sequenced to ensure that the pre-mir sequences match reference sequence imiRBase . The
E. coli Selection:	Kanamycin (25 ug/mL)
Mammalian Cell Selection:	Neomycin
Components:	1 x 10 ug vial miRNA plasmid lyophilized in a 2-D bar-coded Matrix tube
Also Contains:	Both mature and minor (star) miRNA can be over expressed from the construct.
RefSeq:	MI0003482
Synonyms:	rno-mir-215
Locus ID:	100314074



[View online »](#)